

Technical Datasheet

Date: 01.06.17

Page 1 of 2

EN ISO 13485 certified

EN ISO certified		
Article - No:	08-602-6043	Revision: 5
Remarks		single use only

Specification

Product	Universal container PP, 125ml, Ø55x72mm, fitted screw cap PE, 5pcs/bag, container transp., cap yellow, grad., writing area, np pcr ready, sterile R, Autocl. up to 121°C (no cap), CE/IVD
Classification	CE/IVD
Dimensions	upper Ø: 48,0mm (+/- 0,2mm) bottom Ø: 55,0mm (+/- 0,2mm) height (without cap): 71,7mm (+/- 0,2mm) height (w/ cap): 72,0mm (+/- 0,2mm)
Colour	container transp., cap yellow
Graduation	grad., writing area
Nominal Capacity	125ml
Material Characteristics	container: polypropylene (PP) cap: polyethylene (PE)
Temperature range	-40°C-75°C
Autoclavability	Autocl. up to 121°C (no cap)
Product Characteristics	biol. inert physiol. harmless no heavy metals acc. to EC/94/62 Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	free of PCR inhibitors free of human DNA free of RNase and DNase free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Packaging	pcs per bag: 5 pcs per carton: 180 pcs per pallet: 4.320
Storage	dry condition: 40%-60% rel. humidity