

Technical Datasheet

Date: 14.02.17

Page 1 of 1



EN ISO 13485 certified

		EN ISO certified
Article - No:	08-602-7054	Revision: 4
Remarks		single use only

Specification

Product	Universal container PP, 125ml, Ø55x72mm, fitted screw cap PE and label, container transp., cap yellow, grad., np pcr ready, sterile R, Autocl. up to 121°C (no cap), CE/IVD
Classification	CE/IVD
Dimensions	upper Ø: 48,0mm (+/- 0,2mm) bottom Ø: 55,0mm (+/- 0,2mm) height (without cap): 71,7mm (+/- 0,2mm) height (w/ cap): 72,0mm (+/- 0,2mm)
Colour	container transp., cap yellow
Graduation	grad.
Nominal Capacity	125ml
Material Characteristics	container: polypropylene (PP) cap: polyethylene (PE)
Temperature range	-40°C-75°C
Autoclavability	Autocl. up to 121°C (no cap)
Product Characteristics	biol. inert physiol. harmless no heavy metals acc. to EC/94/62 Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	free of PCR inhibitors free of human DNA free of RNase and DNase free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Packaging	pcs per bag: 240 pcs per carton: 240 pcs per pallet: 5.760
Minimum order qty.	240 pcs
Storage	dry condition: 40%-60% rel. humidity bags & cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable