

Technical Datasheet

Date: 04.01.18

Page 1 of 2

| EN ISO certified | | |
|----------------------|--------------------|------------------------|
| Article - No: | 08-602-7802 | Revision: 3 |
| Remarks | | single use only |

Specification

| | |
|---------------------------------|--|
| Product | Universal container PP, 125ml, Ø55x72mm, fitted screw cap PE, seal of originality, container transp., cap blue, grad., writ. area, np pcr ready, sterile R, CE/IVD |
| Classification | CE/IVD |
| Dimensions | upper Ø: 48,0mm (+/- 0,2mm) |
| | bottom Ø: 55,0mm (+/- 0,2mm) |
| | height (w/o cap): 71,7mm (+/- 0,2mm) |
| | height (w/cap): 72,0mm (+/- 0,2mm) |
| Colour | container transp., cap blue |
| Graduation | grad., writ. area |
| Nominal Capacity | 125ml |
| Material Characteristics | cont.: polypropylene (PP) |
| | cap: polyethylene (PE) |
| Temperature range | -40°C-75°C |
| Product Characteristics | biol. inert |
| | physiol. harmless |
| | no heavy metals acc. to EC/94/62 |
| | Food contact: raw material in compl. with reg. EC 1935/2004 |
| Purity Grade | free of PCR inhibitors |
| | free of human DNA |
| | free of RNase and DNase |
| | free of pyrogens |
| Sterilisation Process | radiation sterilization SAL 10-6 |
| Packaging | pcs per bag: 250 |
| | pcs per carton: 250 |
| | pcs per pallet: 5.000 |
| Storage | dry condition: 40%-60% rel. humidity |
| | bags&cases orig. closed |