

Technical Datasheet

Date: 14.02.17

Page 1 of 1

EN ISO certified		
Article - No:	08-602-8051	Revision: 11
Remarks		single use only

Specification

Product	Universal container PP, 125ml, Ø55x72mm, fitted screw cap PE, indiv. packed, container transp., cap red, grad., writing area, np pcr ready, sterile R, Autocl. up to 121°C (no cap), CE/IVD
Classification	CE/IVD
Dimensions	upper Ø: 48,0mm (+/- 0,2mm) bottom Ø: 55,0mm (+/- 0,2mm) height (without cap): 71,7mm (+/- 0,2mm) height (w/ cap): 72,0mm (+/- 0,2mm)
Colour	container transp., cap red
Graduation	grad., writing area
Nominal Capacity	125ml
Material Characteristics	container: polypropylene (PP) cap: polyethylene (PE)
Temperature range	-40°C-75°C
Autoclavability	Autocl. up to 121°C (no cap)
Product Characteristics	biol. inert physiol. harmless no heavy metals acc. to EC/94/62 Food contact: raw material in compl. with reg. EC 1935/2004
Sterilisation Process	radiation sterilization SAL 10-6
Packaging	pcs per cs: 200 pcs per carton: 200 pcs per pallet: 4.800
Minimum order qty.	200pcs
Storage	dry condition: 40%-60% rel. humidity bags & cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable

Approved by nerbe plus quality department