

Technical Datasheet

Date: 11.11.16

Page 1 of 1

		EN ISO certified
Article - No:	08-633-0025	Revision: 1
Remarks		single use only

Specification

Product	Lid PE, without spout, suitable for universal container 08-622-0000, blue
Colour	blue
Material Characteristics	polyethylene (PE)
Temperature range	-50°C - 60°C
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Packaging	pcs per bag: 500 pcs per carton: 1.000 pcs per pallet: 30.000
Minimum order qty.	1.000 pcs
Storage	in dry condition: 40% - 60% rel. humidity bags and cartons originally closed storage temperature +10°C to maximal +30°C pallets cannot be stacked
Approved by nerbe plus quality department	