

# Technical Datasheet

Date: 17.05.17

Page 1 of 1

EN ISO certified		
<b>Article - No:</b>	<b>04-043-0600</b>	<b>Revision: 2</b>
<b>Remarks</b>	<b>fitting tubes: 04-082-0350, 04-082-0550, 04-082-0650, 04-082-9050, 04-083-0096 04-083-0150, 04-083-0640</b>	<b>single use only</b>

## Specification

<b>Product</b>	Cap for PCR microcentrifuge tubes 0,1 ml & 0,2 ml PP, flat, 12-strips, highly transparent, np pcr ready, max. RCF 20.000g, autocl. up to 121°C, CE/IVD
<b>Classification</b>	CE/IVD
<b>Colour</b>	highly transparent
<b>Material Characteristics</b>	polypropylene (PP)
<b>Temperature range</b>	-80°C - 121°C
<b>Autoclavability</b>	autocl. up to 121°C
<b>Centrifugation</b>	max. RCF 20.000g
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
<b>Purity Grade</b>	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
	free of endotoxins
<b>Shelf - Life</b>	not limited at controlled storage conditions
<b>Packaging</b>	pcs per bag: 80
	pcs per carton: 800
	pcs per pallet: 54.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags & cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	