


PDS No. 622201	PRODUCT DATA SHEET				Page 1 of 1
Revision 01	Reaction Tube, 5 ml, Snap Cap				 greiner bio-one
	Greiner Item-No. 622201				
Valid for Item-No.:	622201				

1.	Description / Specification	
1.1	Description	Reaction tube, conical, 5 ml, with attached cap, graduation and writing area
1.2	Dimensions	See Customer Drawing
1.3	Volume	Total volume: 5 ml
1.4	Material / Resin	PP (Polypropylene), free of heavy metal
1.5	Colour	natural
1.6	Sterilization	No
1.7	Quality Control	Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens
2.2	Temperature range	-80°C to +121°C
2.3	Autoclavability	yes
2.4	Centrifugation, max. RCF	25.000 x g: fixed-angle rotor with appropriate support
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	100
3.2	Pieces / Box	1000
3.3	Lot-No.	5-digit 12345 (digit 1 & 2 identify the year and 3 & 4 identify the week)
3.4	Other Information	-

4.	Other Information	
	-	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision	Date 26 March 2018	Date 05 April 2018	Date 09 April 2018	
Date	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.