

Technical Datasheet

Date: 11.11.16

Page 1 of 1

EN ISO certified		
Article - No:	09-302-0020	Revision: 2
Remarks	max. vol.: approx. 3 l / working vol.: approx. 1,25 l	single use only

Specification

Product	Autoclavable bag CPP, 200x300 mm, thickness: 50µm natural, autoclavable up to 134°C
Dimensions	width: 200,00 mm (+/- 6 mm)
	length: 300,00 mm (+/- 15 mm)
	thickness: 50 µm (+/- 10%)
Colour	natural
Material Characteristics	cast polypropylene (CPP)
Autoclavability	autoclavable up to 134°C
Product Characteristics	biologically inert
	physiologically harmless
	side seam
	density approx. 0,90 gr.
Packaging	pcs per carton: 1.000
	pcs per pallet: 60.000
	pcs per dsp case: 100
Minimum order qty.	1.000 pcs
Storage	in dry condition: 40% - 60% rel. humidity
	bags and cartons originally closed
	storage temperature +10°C to maximal +30°C
	pallets cannot be stacked
Approved by nerbe plus quality department	