

Technical Datasheet

Date: 18.05.17

Page 1 of 1

EN ISO certified		
Article - No:	09-322-0000	Revision: 2
Remarks	max. vol.: approx. 30 l / working vol.: approx. 18 l	single use only

Specification

Product	Autoclavable bag PP, 400x800 mm, thickness: 50µm, natural, autoclavable up to 134°C
Dimensions	width: 400,00 mm (+/- 8 mm)
	length: 800,00 mm (-10/+20 mm)
	thickness: 50 µm (+/- 10%)
Colour	natural
Material Characteristics	polypropylene (PP)
Autoclavability	autocl. up to 134°C
Product Characteristics	biologically inert
	physiologically harmless
	extremely resilient double bottom seam
	density approx. 0,90 gr.
Packaging	pcs per cs: 250
	pcs per carton: 250
	pcs per pallet: 15.000
Storage	dry condition: 40%-60% rel. humidity
	bags & cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

Approved by nerbe plus quality department