

Product Specification

RED PEARL

Disposable Nitrile Examination Gloves

- Powder free
- Red
- Chemical resistance protective gloves
- Latex-free
- Suitable for handling food
- Non-sterile



SAFETY MEETS QUALITY

CE 0197



Material

100% Synthetic Nitrile Latex

Performance requirements

Characteristics	Inspection Level	Acceptable Quality Level
Dimensions	S-2	4.0
Physical Properties	S-2	4.0
Freedom from Holes	G-1	1.5
Major Visual defects	G-1	1.5
Minor Visual defects	G-1	4.0

Dimensions

Number	9801	9802	9803	9804	9805
Size	XS = 5–6	S = 6–7	M = 7–8	L = 8–9	XL = 9–10
Width	75 ± 5 mm	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	≥ 110 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Thickness

Measuring Points	Single Wall (mm)
Cuff (distance 25 mm to the beaded edge)	0.08 ± 0.03
Palm (in the middle of the palm)	0.10 ± 0.03
Fingertip (distance 15 mm to top of middle finger)	0.12 ± 0.03

Physical Properties

Standard	Before Aging	After Aging
Tensile Strength (N)	min. 6.0	min. 6.0

Powder

No powder lubricant added

Protein content

No proteins contained

Cuff

Beaded cuff

Design

Ambidextrous

Surface

Micro-roughened finger tips

Color

Red

Quality assurance

EG-construction sample test TÜV Rheinland – CE 0197

PSA – EN 374 CAT III, EN 388, EN 420, 89/686/EEC

MPG – EN 455 Teil 1–4

EG-Standard 596/2009

EG-Standard 202/2014

EN ISO 9001:2008

EN ISO 13485:2012

Packing Style

100 pcs per dispenser, 10 dispensers per carton

Hint

The above information apply to the standard version.
Subject to change.

HUBERLAB. AG
Industriestrasse 123
4147 Aesch

T 061 717 99 77
F 061 711 93 42

www.huberlab.ch
info@huberlab.ch

HUBERLAB.

committed to science