

Product Specification

CLASSIC

Disposable Latex Examination Gloves

- Prepowdered
- White
- Non-sterile



SAFETY MEETS QUALITY

Material

Natural High Grade Rubber Latex

Performance requirements

Characteristics	Inspection Level	Acceptable Quality Level
Dimensions	S-2	4.0
Physical Properties	S-2	4.0
Freedom from Holes	G-1	1.5
Major Visual defects	G-1	2.5
Minor Visual defects	G-1	4.0

Dimensions

Number	1101	1102	1103	1104	1105
Size	XS = 5-6	S = 6-7	M = 7-8	L = 8-9	XL = 9-10
Width	75 ± 5 mm	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	≥ 110 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Thickness

Measuring Points	Single Wall (mm)
Cuff (distance 25 mm to the beaded edge)	0.11 ± 0.03
Palm (in the middle of the palm)	0.13 ± 0.03
Fingertip (distance 15 mm to top of middle finger)	0.14 ± 0.03

Physical Properties

Standard	Before Aging	After Aging
Median force at break (N)	min. 6.0	min. 6.0
Ultimate Elongation (%)	min. 650	min. 500
Tensile Strength (MPa)	min. 18.0	min. 14.0

Powder

Absorbable cornstarch according USP

Cuff

Beaded cuff

Design

Ambidextrous

Surface

Smooth

Color

White

Quality assurance

Comply the requirements of ASTM D3577-91, EN 455 part 1-2 + 4, Manufactured under ISO 9001:2008, ISO 13485:2012, according GMP

Packing Style

100 pcs per dispenser, 10 dispensers per carton

Hint

The above information apply to the standard version.
Subject to change.