

Plast 4 Labs



HUBERLAB.

committed to science

Index

Lab Bottles

WIDEMOUTH BOTTLES RECTANGULAR BASED - GRADUATED	Pg 2
WIDE MOUTH BOTTLES - GRADUATED	Pg 2
NARROW MOUTH BOTTLES - GRADUATED	Pg 3
NARROW MOUTH BOTTLES - GREY	Pg 3
PP WIDE NECK BOTTLES - GL THREAD	Pg 4
PE WIDE NECK BOTTLES - GL THREAD	Pg 4
PP WIDE NECK BOTTLES WITH GRADUATION - GL THREAD	Pg 5
PE WIDE NECK BOTTLES AMBER - GL THREAD	Pg 5
PP NARROW NECK BOTTLES - GL THREAD	Pg 6
PE NARROW NECK BOTTLES AMBER - GL THREAD	Pg 6
PE NARROW NECK BOTTLES - GL THREAD	Pg 7
GL CAPS WITH TUBES	Pg 7
GL CAPS	Pg 7
GL CAPS LONG DROPPER	Pg 7
GL CAPS SHORT DROPPER	Pg 7
PP WIDE NECK BOTTLES	Pg 8
PP NARROW NECK BOTTLES	Pg 8
HEXAGONAL CONIC CAPS	Pg 8
ROUND BOTTLE WITH SEALING CAP	Pg 9
APPROVED UN BOTTLES	Pg 9
NARROW MOUTH BOTTLES WITH SEALING CAP	Pg 9
WIDE MOUTH BOTTLES WITH SEALING CAP	Pg 10
BOTTLES WITH CHILD PROOF CAP	Pg 10
BOTTLES WITH DISPENSER	Pg 10
NARROW MOUTH BOTTLES SQUARE BASED	Pg 11
WIDEMOUTH BOTTLES SQUARE BASED	Pg 11
NARROW MOUTH BOTTLES - SEALING CAP	Pg 11
WIDE MOUTH BOTTLES - SHORT TYPE	Pg 12
WIDE MOUTH BOTTLES SCREW CAP	Pg 12
SPRAY BOTTLES	Pg 12
DROP BOTTLES	Pg 13
WASH BOTTLES	Pg 13
ADHESIVE LABELS	Pg 13
WASH BOTTLES WITH INTEGRATED TUBE	Pg 14
WIDE MOUTH BOTTLES SCREW CAP	Pg 14
WIDE MOUTH BOTTLE SEALING CAP	Pg 14

Jars

JARS WITH SCREW CAP AND INNER CAP	Pg 16
JARS WITH SCREW CAP AND INNER CAP	Pg 16
JARS WITH PRESSURE SEALING CAP	Pg 16
SQUARE BASED JARS	Pg 17
STACKABLE SQUARE BASED JARS	Pg 17
SQUARE BASED JAR	Pg 17
CYLINDRICAL JARS - PRESSURE CAP	Pg 18
HERMETIC JARS PRESSURE CAP	Pg 18
HERMETIC JARS SCREW CAP	Pg 18
CONICAL JARS	Pg 19
TRANSPARENT JARS	Pg 19
BOTTLE WITH BRUSH	Pg 19

Carboy - Aspirators - Storage Bottles

CARBOY/ ASPIRATOR	Pg 22
WIDE MOUTH DRUMS	Pg 22
STORAGE BOTTLES	Pg 23
TANKS	Pg 24
TANKS WITH STOPCOCK	Pg 24
UN APPROVED TANKS	Pg 25
WIDE MOUTH DRUM	Pg 25
WIDE MOUTH BOTTLE WITH HANDLE	Pg 25

Funnels - Vials - Graduated Jugs

FUNNELS FOR TRANSFER	Pg 28
FUNNEL WITH FILTER	Pg 28
FUNNELS FOR POWDER	Pg 28
FUNNELS FOR ANALYSIS	Pg 29
FUNNELS WITH RIBS	Pg 29
SCINTILLATION VIALS	Pg 29
JUGS SHORT TYPE	Pg 30
JUGS TALL TYPE	Pg 30
JUGS TALL TYPE WITH BLUE GRADUATION	Pg 30

Cylinders - Beakers - Flasks

TALL TYPE CYLINDERS	Pg 32
TALL TYPE CYLINDERS WITH BLUE GRADUATION	Pg 32
SHORT TYPE CYLINDERS	Pg 33
BEAKERS	Pg 33
CONICAL FLASKS ERLLENMEYER	Pg 34
VOLUMETRIC FLASKS	Pg 34
CONICAL CAPS	Pg 34

Tubing - Monobulb - Ballons - Stoppers - Stirring bars

PVC TUBING	Pg 36
SYLICON TUBING	Pg 36
PTFE TUBING	Pg 36
RED RUBBER TUBING	Pg 36
LATEX TUBING	Pg 36
PARA TUBING	Pg 36
MONO BULB	Pg 37
VACUUM BULB	Pg 37
RICHARDSON BULBS	Pg 37
GAS SAMPLING BALLOONS	Pg 37
REVERSIBLE STOPPERS	Pg 37
SYLICON CONICAL STOPPERS	Pg 37
RED RUBBER STOPPERS	Pg 37
RED RUBBER STOPPERS WITH HOLES	Pg 37
MAGNETIC STIRRING BARS	Pg 37
LATEX BULB	Pg 37

Containers - Scoops

STACKABLE CRATES	Pg 40
ISOTHERMAL CONTAINERS	Pg 42
SCOOPS	Pg 42
LARGE CUP	Pg 42
SPATULAS	Pg 42
PUMP	Pg 43
BOTTLES CARRIER	Pg 44
BED BASIN	Pg 44
MALE BED BOTTLE	Pg 44

Lab bottles




Bottles

POLYETHYLENE BOTTLES SUITABLE FOR THE STORAGE OF SAMPLES: LIQUIDS, POWDERS AND SOLIDS. THEY ARE EQUIPPED WITH CAP AND INNER-CAP.

On the shoulder and cap there are holes, that allow the use of seals or slips.

The depth of the thread of the neck and the cap, together with the presence of the inner-cap, ensure a perfect seal. The translucence and graduation allow a good view of the content, while the complete emptying is guaranteed by the inclination of the shoulders.

They are supplied with inner-cap and cap and packed in cardboard boxes.



WIDE MOUTH BOTTLES RECTANGULAR BASED - GRADUATED

Material: High density Polyethylene

Color: Neutral – For food use

The rectangular shape allows a better use of the spaces.

The precise graduation allows the filling without the use of beakers or cylinders. They are equipped with inner-cap to ensure a perfect seal. Suitable for the storage of dense fluids, powders, granulates.

Confidential notes can be placed in the space between cap and inner-cap.

Item	Capacity ml	Graduation ml	Dimension mm	Height mm	Mouth mm
15.1608.06	25	5	30x32	53	18,5
15.1609.01	50	10	35x37	70	24
15.1610.08	100	20	44x47	90	34,5
15.1611.01	250	50	57x60	110	34,5
15.1612.06	500	100	70x80	138	45
15.1613.01	1000	100	82x102	176	60
15.1614.02	2000	100	100x134	210	60



WIDE MOUTH BOTTLES - GRADUATED

Material: Low density Polyethylene.

Color: Neutral – For food use

They are equipped with inner-cap for a better seal and with holes on the shoulder and on the cap to enable the use of security seals.

The precise graduation allows the filling without the use of beakers or cylinders. The total emptying is ensured by the inclination of the shoulders of the bottle. Suitable for the storage of dense fluids, powders, granulates. Confidential notes can be placed in the space between cap and inner-cap.

Item	Capacity ml	Graduation ml	Diameter mm	Height mm	Mouth mm
15.1408.01	50	10	40	77	24
15.1409.01	100	20	48	89	24
15.1410.03	250	25	61	126	34,5
15.1411.01	500	50	75	161	34,5
15.1412.01	1000	100	93	201	45
15.1413.01	2000	100	115	250	45

POLYETHYLENE BOTTLES SUITABLE FOR THE STORAGE OF SAMPLES: LIQUIDS, POWDERS AND SOLIDS. THEY ARE EQUIPPED WITH CAP AND INNER-CAP.

On the shoulder and cap there are holes, that allow the use of seals or slips.

The depth of the thread of the neck and the cap, together with the presence of the inner-cap, ensure a perfect seal. The translucence and graduation allow a good view of the content, while the complete emptying is guaranteed by the inclination of the shoulders.

They are supplied with inner-cap and cap and packed in cardboard boxes.

NARROWMOUTH BOTTLES - GRADUATED

Material: Low density Polyethylene

Color: Neutral - For food use

They are equipped with inner-cap for a perfect seal and with holes on the shoulder and on the cap to allow the use of seals.

Suitable for the storage of very fluid liquids. The graduation precise allows the filling without the use of beakers or cylindres.

Item	Capacity ml	Graduation ml	Diameter mm	Height mm	Mouth mm
15.1317.01	50	10	39	76	18,5
15.1319.01	100	20	48	115	18,5
15.1323.02	250	25	60	138	23
15.1324.02	500	50	74	165	23
15.1325.02	1000	100	93	213	34,5
15.1326.02	2000	100	116	273	34,5



NARROWMOUTH BOTTLES - GREY

Material: Low density Polyethylene

Color: Grey - For food use

Essential for photosensitive products, are equipped with inner-cap for a better seal and with holes on the shoulder and on the cap to allow the use of seals.



Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1323.03	250	60	138	23
15.1324.03	500	74	165	23
15.1325.03	1000	93	213	34,5
15.1326.03	2000	116	273	34,5

WARRANTY SEAL FOR BOTTLES PL 400
Stainless steel cable, 7 wire length 150 mm Progressively numbered. Usable on all bottles with a cap with a hole for seal: Narrow mouth - Wide mouth - Rectangular.



Bottles

General purpose bottles, wide mouth. long neck thread ensures leakproof closure, wide neck for an easy filling and emptying of liquids and powder. Compliant to DIN 13316 and 168. Perfect for samples collections, storage and transport, ideal for liquids, pastes and granular substances. The translucent and graduation allow a good view of the content, while the complete emptying is guaranteed by the inclination of the shoulders.



PP WIDE NECK BOTTLES - GL THREAD

Material: Polypropylene
Color: Neutral - for food use
Can be filled with hot substances
Autoclavable up to 121° C.

Item	Capacity ml	Diameter mm	Neck diameter mm	Height mm	Cap size GL
15.9126.01	50	40	24	77	32
15.9326.01	100	48	24	92	32
15.9526.01	250	64	31	123	40
15.9726.01	500	75	38,5	153	50
15.9926.01	1000	94	50	202	65
15.9136.01	2000	120	50	243	65



PE WIDE NECK BOTTLES - GL THREAD

Material: Low density Polyethylene
Color: Neutral - For food use

Item	Capacity ml	Diameter mm	Neck diameter mm	Height mm	Cap size GL
15.9806.01	50	40	24	77	32
15.9016.01	100	48	24	92	32
15.9216.01	250	64	31	123	40
15.9416.01	500	75	38,5	153	50
15.9616.01	1000	94	50	202	65
15.9816.01	2000	120	50	243	65

General purpose bottles, wide mouth, long neck thread ensures leakproof closure, wide neck for an easy filling and emptying of liquids and powder. Compliant to DIN 13316 and 168. Perfect for samples collections, storage and transport, ideal for liquids, pastes and granular substances. The translucent and graduation allow a good view of the content, while the complete emptying is guaranteed by the inclination of the shoulders.

PP WIDE NECK BOTTLES WITH GRADUATION - GL THREAD

Material: Polypropylene
Color: Neutral - For food use
Can be filled with hot substances
Autoclavable up to 121° C.
With graduation easy to read

Item	Capacity ml	Diameter mm	Neck diameter mm	Height mm	Capsize GL
15.3608.01	50	40	24	77	32
15.3610.01	100	48	24	92	32
15.3612.01	250	64	31	123	40
15.3614.01	500	75	38,5	153	50
15.3616.01	1000	94	50	202	65
15.3618.01	2000	120	50	243	65



PE WIDE NECK BOTTLES AMBER - GL THREAD

Material: Low density Polyethylene
Color: Amber - For food use
Light protection for sensitive substances

Item	Capacity ml	Diameter mm	Neck diameter mm	Height mm	Capsize GL
15.2608.01	50	40	24	77	32
15.2610.01	100	48	24	92	32
15.2612.01	250	64	31	123	40
15.2614.01	500	75	38,5	153	50
15.2616.01	1000	94	50	202	65



Bottles

Bottles for liquids very fluid with narrow neck. There are no shoulders to restrict the flow. Compliant to DIN 13316 and 168. Perfect for sample collection, storage and transport. The translucent and graduation allow a good view of the content, while the complete emptying is guaranteed by the inclination of the shoulders.



PP NARROW NECK BOTTLES - GL THREAD

Material: Polypropylene
Color: Neutral - For food use
Can be filled with hot substances
Autoclavable up to 121° C.

Item	Capacity ml	Diameter mm	Neck diameter mm	Height mm	Cap size GL
15.9595.01	20	32	9	57	14
15.9795.01	100	45	13	102	18
15.9995.01	250	60	18	138	25
15.9106.01	500	75	13	169	25
15.9306.01	1000	95	21	212	28



PE NARROW NECK BOTTLES AMBER - GL THREAD

Material: Low density Polyethylene
Color: Amber - For food use
Light protection for sensitive substances

Item	Capacity ml	Diameter mm	Neck diameter mm	Height mm	Cap size GL
15.2620.01	100	45	13	102	18
15.2622.01	250	60	18	138	25
15.2624.01	500	75	18	169	25
15.2626.01	1000	95	21	212	28

Bottles

Bottles for liquids very fluid with narrow neck. There are no sholders to restrict the flow. Compliant to DIN 13316 and 168
Perfect for sample collection, storage and transport. The translucent and graduation allow a good view of the content, while the complete emptying is guaranteed by the inclination of the shoulders.

PENARROWNECK BOTTLES - GL THREAD

SUPPLIED WITHOUT CAPS

Material: Low density Polyethylene.

Color: Neutral - For food use.

Multipurpose bottles for different use.

They can be sampling bottles, wash bottles, drop bottles.

With long or short dropper.

Item	Capacity ml	Diameter mm	Neck diameter mm	Height mm	Capsize GL
15.9108.01	50	40	13	84	18
15.9167.95	100	45	13	102	18
15.9167.99	250	60	18	138	25
15.9168.01	500	75	18	169	25
15.9188.01	1000	95	21	212	28
15.9109.02	2000	120	21	260	28



GL CAPS WITH TUBES

Item	GL
15.9003.01	18
15.9103.01	25
15.9203.01	28

GL CAPS

Material: High density Polyethylene.

Item	GL
15.9008.01	18
15.9048.01	25
15.9088.01	28



GL CAPS LONG DROPPER

Item	GL
15.9603.01	18
15.9703.01	25
15.9803.01	28

GL CAPS SHORT DROPPER

Item	GL
15.9303.01	18
15.9403.01	25



Bottles

Bottles with pressure cap with wide or narrow neck, the complete emptying is guaranteed by the inclination of the shoulders. Can be filled with hot substances and can be sterilized in autoclave up to 121°C. The translucent and graduation allow a good view of the content, while the complete emptying is guaranteed by the inclination of the shoulders.



PP WIDE NECK BOTTLES

Supplied with caps
Material: Polypropylene
Color: Neutral - For food use


Item	Capacity ml	Diameter mm	Neck diameter mm	Height mm	Cap NS
15.3710.01	100	52	29	96	29/25
15.3712.01	250	72	34	132	34/35
15.3714.01	500	86	45,5	162	45/40
15.3716.01	1000	107	60	200	60/46



PP NARROW NECK BOTTLES

Supplied with caps
Material: Polypropylene
Color: Neutral - For food use

Item	Capacity ml	Diameter mm	Neck diameter mm	Height mm	Cap NS
15.3720.01	100	52	14,5	100	14/23
15.3722.01	250	70	19	128	19/26
15.3724.01	500	86	24	162	24/29
15.3726.01	1000	106	29	197	29/32



HEXAGONAL CONICAL CAPS

Item	Size
15.2240.07	7/16
15.2241.03	10/19
15.2242.01	12/21
15.2243.01	14/23
15.2244.01	14/15

Item	Size
15.2245.01	18/26
15.2246.01	19/26
15.2247.01	24/29
15.2248.01	24/20
15.2249.01	29/32

Item	Size
15.2250.01	29/22
15.2251.01	34/35
15.2252.01	45/40
15.2253.01	60/46
15.2254.01	85/65

Bottles with self sealing cap

ROUND BOTTLE WITH SEALING CAP Ø 41 MOUTH (DIN 50)

Material: High density Polyethylene
Color: White milk - For food use
Self sealing cap

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1435.01	1000	90	220	41



APPROVED U.N. BOTTLES FOR TRANSPORTATION OF DANGEROUS GOODS

Material: High density Polyethylene
Color: White milk - For food use
U.N. omologation 1H1/X1,3 - Self sealing cap
Mouth Ø41mm (DIN 50)

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1451.02	250	63	127	41
15.1452.01	500	69	186	41
15.1453.02	1000	89	240	41



NARROW MOUTH BOTTLES WITH SEALING CAP

Material: High density Polyethylene
Color: Neutral - For food use
Screw cap with opening-off seal. its use is recommended when you need to preserve the content after the filling. They are suitable for the storage of very fluid liquids. The complete emptying is guaranteed by to the inclination of the shoulders.

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1351.02	50	36	84	17,5
15.1352.02	100	44	103	17,5
15.1353.02	250	60	136	22
15.1354.03	500	74	170	22
15.1356.01	1000	90	222	22



Bottles with self sealing cap



WIDE MOUTH Ø 30 MM BOTTLES WITH SEALING CAP

Material: High density Polyethylene.

Color: Neutral - For food use.

Screw cap with opening-off seal. Its use is recommended when you need to preserve the contents after filling. Suitable for the storage of dense fluids, granulates and powders. The complete emptying is guaranteed by the inclination of the shoulders.

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1343.02	250	60	140	30
15.1344.01	500	74	179	30
15.1345.02	750	81	214	30
15.1346.01	1000	90	218	30



BOTTLES WITH SELF SEALING CHILD PROOF CAP

Material: High density Polyethylene

Color: Neutral - For food use.

Cap can be unscrewed only by pushing and turning. The complete emptying is guaranteed by the inclination of the shoulders. They are suitable for the storage of very fluid liquids.

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1373.01	250	60	140	22
15.1374.01	500	74	179	22
15.1376.01	1000	90	218	22



BOTTLES WITH DISPENSER

Material: High density Polyethylene

Color: Neutral - For food use

They are equipped with a no self-sealing dispenser cap and little cap to close the tip.

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1363.01	250	60	136	22
15.1364.01	500	74	170	22
15.1366.01	1000	90	222	22

Bottles with self sealing cap

NARROW MOUTH GRADUATED BOTTLES SQUARE BASE –SEALING CAP

Material: High density Polyethylene

Color: Neutral - For food use

They are equipped with inner-cap. The squared shape allows a better use of the spaces. Graduation is in embossed form, which allows the filling without the use of beakers or cylinders.

Suitable for the storage of very fluid liquids.

Item	Capacity ml	Graduation ml	Dimension mm	Height mm	Mouth mm
15.1523.48	250	50	50x50	140	21
15.1525.01	500	100	67x67	165	21
15.1526.01	1000	100	80x80	210	21



WIDE MOUTH GRADUATED BOTTLES SQUARE BASE SEALING CAP

Material: High density Polyethylene

Color: Neutral - For food use

The squared shape allows a better use of the spaces. Graduation is in embossed form, which allows the filling without the use of beakers or cylinders. Suitable for the storage of dense fluids, granulates and powders.

Item	Capacity ml	Graduation ml	Dimension mm	Height mm	Mouth mm
15.1524.01	250	50	50x50	140	30
15.1535.01	500	100	65x65	165	30
15.1536.01	1000	100	80x80	210	30



NARROW MOUTH BOTTLES – SEALING CAP

Material: High density Polyethylene

Color: Neutral - For food use.

They are equipped with inner-cap. They are also available with perforated inner-cap: minimum 1000 PCs. The self-sealing cap is tamper proof of the substance introduced. Suitable for the storage of very fluid liquids

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
• 15.1380.02	25	27	69	12
• 15.1381.02	50	31	88	12
• 15.1382.02	100	35	117	12
15.1383.02	250	52	167	20,5
15.1384.04	500	64	200	20,5
15.1385.02	1000	86,5	220	20,5

• Not self seal cap



Bottles



WIDE MOUTH BOTTLES – SHORT TYPE

Material: High density Polyethylene

Color: Neutral - For food use

Cap with gasket. Emptying guaranteed by the inclination of the shoulders. Thanks to their reduced height they can be placed in drawers or very low spaces. Suitable for the storage of very fluid liquids.

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1302.03	100	50	97	29
15.1303.03	150	60	97	29
15.1304.03	200	70	97	29



WIDE MOUTH BOTTLES WITH SCREW CAP

Material: High density Polyethylene

Color: Neutral - For food use

Equipped with inner-cap. Emptying guaranteed by the inclination of the shoulders. The wide mouth allows an easy filling. Suitable for the storage of dense liquids, granulates and powders.

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1035.01	1300	100	215	60
15.1036.01	2000	120	245	60
15.1037.03	3000	140	280	75

- Cap with Polespan gasket



SPRAY BOTTLES

Material: Low density Polyethylene

Color: Neutral - For food use

Item	Capacity ml
15.1394.03	500
15.1395.02	900
15.1396.03	spray parts

Drop and wash bottles

DROP BOTTLES

Material: Low density Polyethylene

Color: Neutral - For food use

Item	Capacity ml	Diameter mm	Height mm
15.1101.02	50	35	118
15.1102.03	125	44	188
15.1103.03	250	58	220
15.1105.01	500	70	253
15.1106.03	750	81	280
15.1107.06	1000	90	295



WASH BOTTLES

Material: Low density Polyethylene

Color: Neutral - For food use

They have excellent clarity and easy use through a slight pressure on the body of the bottle. Tube with removable dispensing tip.

Item	Capacity ml	Diameter mm	Height mm
15.1180.01	100	42	218
15.1182.01	250	58	218
15.1183.01	500	70	218
15.1184.01	750	81	253
15.1185.01	1000	90	270



ADHESIVE LABELS FOR WASH BOTTLES

Dimension mm 130x35

Item	Description	Color	Pack PCs
15.8900.01	Dist. Water	Blue	10
15.8900.02	Acetone	Red	10
15.8900.03	Methanol	Black	10
15.8900.04	Isopropanol	Brown	10
15.8900.05	Ethyl Acetate	Orange	10
15.8900.06	Ethanol	Green	10
15.8900.07	Chloroform	Black	10
15.8900.08	N-Hexane	Black	10
15.8900.09	Methylene Chloride	Black	10
15.8900.10	Toluol	Black	10
15.8901.01	Xylol	Black	10



Wash bottles and bottles



WASH BOTTLES WITH INTEGRAL TUBE

Material: Low density Polyethylene

Color: Neutral - Capacity: 500 ml

Wash bottles with wide mouth for easy filling. Dispensing tube fully molded with the cap. The name of the solution inside the bottle is printed in: German, English, French and Spanish.

Item	Solution name
15.1470.08	acetone
15.1471.01	distilled water
15.1472.08	ethanol
15.1473.01	ethyl acetate
15.1474.01	isopropanol
15.1475.01	methanol
15.1476.01	methylethylketone
15.1477.01	tetrahydrofuran
15.1478.01	toluene



WIDE MOUTH BOTTLES WITH SCREW CAP

Material: High density Polyethylene

Color: Neutral - For food use

Equipped with inner cap and cap with holes for warranty seal.

Emptying guaranteed by the inclination of the shoulders.

The wide mouth allows an easy filling. Suitable for the storage of dense liquids, granulates and powders.

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1031.01	2000	125	235	82
15.1032.03	3000	135	285	82



WIDE MOUTH BOTTLE WITH SEALING CAP

Material: High density Polyethylene

Color: Neutral - For food use

Self-sealing cap DIN51 with gasket. The self-sealing cap is tamper proof of the substance introduced.

Minimum for the production 5000 PCs.

Item	Capacity ml	Graduation ml	Diameter mm	Height mm	Mouth mm
15.1448.01	2000	100	115	250	45

Jars



Jars with screw and pressure cap



JARS WITH SCREW CAP AND INNER-CAP

Material: High density Polyethylene

Color: Neutral - For food use

Thanks to the wide mouth they are suitable for storage: creams, semi-solid substances, crystals, granulates and powders.

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1020.05	70	50	60	36
15.1021.01	120	57	70	36
15.1022.01	250	67	90	50



JARS WITH SCREW CAP AND INNER-CAP

Material: High density Polyethylene

Color: Neutral - For food use

Thanks to the wide mouth they are suitable for storage: creams, semi-solid substances, crystals, granulates and powders.

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1023.02	500	91	103	70
15.1025.10	750	96	120	70
15.1024.01	1000	110	127	86
15.1026.02	2000	125	185	86



JARS WITH PRESSURE SEALING CAP

Material: High density Polyethylene

Color: Neutral - For food use

The self-sealing pressure cap is tamper proof. Thanks to the wide mouth they are suitable for the storage of creams, semi-solid substances, granulates and powders.

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1054.01	500	105	74	86
15.1027.01	1000	110	130	86
15.1028.02	1500	110	180	86
15.1029.01	2000	125	185	86

• Color: White milk

Jars with self-sealing cap

SQUARE BASED JARS WITH PRESSURE SEALING CAP

Material: High density Polyethylene

Color: White milk - For food use

The self-sealing pressure cap is tamper proof. Thanks to the wide mouth they are suitable for the storage of creams, semi-solid substances, granulates and powders. The square shape allows a better use of the spaces.

Item	Capacity ml	Dimension mm	Height mm	Mouth mm
15.1043.02	700	90x90	115	70
15.1045.02	1300	120x120	134	85



STACKABLE SQUARE BASE JARS, WITH SCREW SELF-SEALING CAP

Material: High density Polyethylene

Color: Neutral - For food use

The square shape allows a better use of the spaces.

The self-sealing screwing cap is tamper proof. Thanks to the wide mouth they are suitable for the storage of creams, semi-solid substances, granulates. and powders.

Item	Capacity l	Dimension mm	Height mm	Mouth mm
15.1055.43	5	195x195	215	110
15.1058.02	10	195x195	400	110



SQUARE BASED JAR - SCREW SELF-SEALING CAP

Material: High density Polyethylene

Color: Neutral - For food use

Item	Capacity ml	Dimension mm	Height mm	Mouth mm
15.1539.01	1000	95x95	135	60



Jars with self-sealing cap



CYLINDRICAL JAR SELF-SEALING PRESSURE CAP

Material: Polypropylene

Color: White milk - For food use

The screw cap ensure pratical use for a perfect opening an closing of the jar. Thanks to the wide mouth they are suitable to storage: creams, semi-solid substances, crystals, granulates and powders. "

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.2195.01	650	105	100	95



HERMETIC JARS - SELF-SEALING AND PRESSURE CAP

Material: Polypropylene

Color: White milk - For food use

They guaranteed the tightness of the content once you break the seal. Jars can be closed again after the first time by pushing the cap on the mouth. Thanks to the wide mouth they are suitable to storage: creams, semi-solid substances, crystals, granulates and powders.

Item	Capacity ml	Diameter mm	Height mm	Cap code
15.2415.01	50	35	70	15.2415.02
15.2418.01	100	45	80	15.2418.02
15.2420.01	150	55	70	15.2420.02
15.2426.01	280	65	100	15.2426.02
15.2430.01	500	90	85	15.2430.02
15.2433.01	1000	110	115	15.2433.02
15.2434.01	1500	110	185	15.2434.02

To purchase the complete jar, please ask for body and for cap code



HERMETIC JARS - SELF-SEALING SCREW CAP

Material: Polypropylene

Color: White milk - For food use

The screw cap ensure pratical use for a perfect opening and closing of the jar. Thanks to the wide mouth they are suitable to storage: creams, semi-solid substances, crystals, granulates and powders.

Item	Capacity ml	Diameter mm	Height mm	Cap code
15.2450.01	50	35	50	15.2450.02
15.2452.01	100	50	60	15.2452.02
15.2454.01	250	65	85	15.2454.02
15.2458.01	500	90	110	15.2458.02
15.2462.01	1000	110	120	15.2462.02
15.2464.01	1250	110	155	15.2464.02

To purchase the complete jar, please ask for body and for cap code

CONICAL JARS WITH PRESSURE SELF-SEALING CAP

Material: Polypropylene

Color: White milk - For food use

Self-sealing cap. Thanks to the wide mouth they are suitable to storage: creams, semi-solid substances, crystals, granulates and powders.

Item	Capacity ml	Diameter mm	Height mm
15.1242.01	200	86	55
15.1243.01	300	86	81
15.1245.01	650	113	92
15.1246.01	1000	113	133



TRANSPARENT JARS (CYLINDRICAL)

Material: Crystal Polystyrene

For food use

Transparent body with low density polyethylene pressure cap.

Thanks to the wide mouth they are suitable to contain: creams, semi-solid substances, crystals, granulates and powders.

Item	Capacity ml	Diameter mm	Height mm
15.2174.01	30	39	32
15.2175.05	50	39	52
15.2177.01	100	51	60
15.2178.01	150	65	50
15.2180.01	200	65	67
15.2182.01	250	65	85
15.2184.01	500	100	95

• Conical



GRADUATED BOTTLE WITH BRUSH AND STEEL BALL

Material: High density Polyethylene

Color: Neutral, green cap - For food use

Item	Capacity ml	Diameter mm	Height mm	Mouth mm
15.1307.01	50	40	73	30





HUBERLAB.

committed to science

Carboy-aspirator Storage bottles



Carboy and aspirator



CARBOY/ASPIRATOR

Material: High density Polyethylene - Color: Neutral - For food use
Very strong and durable for extended use over time with carrying handle. On request we can provide stopcock BSP 3/4" threaded with a depth of 12 mm. with gasket for a perfect seal. They can be used with temperatures below 0°C. Cap with holes for warranty seal.



STOPCOCK 15.1393.02
Material: High density Polyethylene 3/4" threaded
For food use.
Outlet hole ø 9 mm.

CARBOY (WITHOUT STOPCOCK)

Item	Capacity lt	Diameter mm	Height mm	Mouth mm
15.1187.01	5	167	355	51
15.1188.01	10	212	430	51

ASPIRATOR (WITH STOPCOCK)

Item	Capacity lt	Diameter mm	Height mm	Mouth mm
15.1187.02	5	167	355	51
15.1188.02	10	212	430	51



CARBOY/ASPIRATOR

Material: High density Polyethylene - Color: Neutral - For food use
Very strong and durable for extended use over time with carrying handle. On request we can provide stopcock BSP 3/4" threaded with a depth of 12 mm. with gasket for a perfect seal. They can be used with temperatures below 0°C. Cap with holes for warranty seal.



STOPCOCK 15.1393.02
Material: High density Polyethylene 3/4" threaded
For food use.
Outlet hole ø 9 mm.

CARBOY (WITHOUT STOPCOCK)

Item	Capacity lt	Diameter mm	Height mm	Mouth mm
15.1190.50	15	240	465	51
15.1190.55	20	265	510	51
15.1191.01	25	275	580	80
15.1192.01	50	350	715	80

ASPIRATOR (WITH STOPCOCK)

Item	Capacity lt	Diameter mm	Height mm	Mouth mm
15.1190.00	15	240	465	51
15.1190.02	20	265	510	51
15.1191.02	25	275	580	80
15.1192.02	50	350	715	80



WIDE MOUTH DRUMS

Material: High density Polyethylene - Color: Neutral - For food use
Very strong and durable for extended use over time with carrying handle. On request we can provide stopcock BSP 3/4" threaded with a depth of 12 mm. with gasket for a perfect seal. They can be used with temperatures below 0°C. Cap with holes for warranty seal.



STOPCOCK 15.1393.02
TO BE ORDERED SEPARATELY
Material: High density Polyethylene 3/4" threaded
For food use.
Outlet hole ø 9 mm.

Item	Capacity lt	Diameter mm	Height mm	Mouth mm
15.1386.02	30	380	440	250
15.1387.02	50	380	650	250
15.1388.02	100	450	800	250

Storage bottles

STORAGE BOTTLES (ECONOMIC TYPE)

Material: High density Polyethylene

Color: Neutral - For food use

They are equipped with inner-cup

They are equipped with inner-cup for a tight seal.

Ideal for containing liquids, powders and granulates.

Item	Capacity lt	Diameter mm	Height mm	Mouth mm
15.1505.01	5	160	280	88
15.1510.06	10	193	350	88



WIDE MOUTH STORAGE BOTTLES

Material: High density Polyethylene

Color: Neutral - For food use

They are equipped with inner-cap and cap with holes for warranty seal.

Very strong and durable for extended use over time. They are equipped with inner-cap for a tight seal. Supplied with cap with holes for warranty seal. Ideal for containing liquids, powders and granulates.

Item	Capacity lt	Diameter mm	Height mm	Mouth mm
15.1414.01	5	167	350	85
15.1415.01	10	210	420	100



WIDE MOUTH STORAGE BOTTLES

Material: High density Polyethylene

Color: Neutral - For food use

They are equipped with inner-cap and cap with holes for warranty seal.

Very strong and durable for extended use over time. They are equipped with inner-cap for a tight seal. Supplied with cap with holes for warranty seal. Ideal for containing liquids, powders and granulates.

Item	Capacity lt	Diameter mm	Height mm	Mouth mm
15.1416.01	15	240	460	115
15.1417.01	20	265	500	130
15.1350.02	25	275	570	130



Tanks



TANKS

Material: High density Polyethylene

Color: Neutral - For food use

Self-sealing cap

Tanks with capacity 1-2 lt are provided with inner-cap

Item	Capacity lt	Dimension mm	Height mm	Mouth mm
15.3941.01	1	69x118	183	21
15.3942.01	2	79x143	235	21
• 15.3943.01	3	95x160	260	21

• Without self-sealing cap



TANKS

Material: High density Polyethylene

Color: Neutral - For food use

Self-sealing cap

Supplied with pouring spout

Item	Capacity lt	Dimension mm	Height mm	Mouth mm
• 15.3948.01	5	120x225	280	44
15.3949.01	10	145x275	340	44
15.3958.01	15	160x315	390	44
15.3950.04	20	170x340	425	44
15.3951.01	25	195x355	500	44
15.3952.01	30	195x355	550	44

• Without pouring spout



TANKS WITH SCREW STOPCOCK

Material: High density Polyethylene

Color: Neutral - For food use

Self-sealing cap - Supplied with Plastirob stopcock **15.5953.01**

Item	Capacity lt	Dimension mm	Height mm	Mouth mm
15.3961.01	10	150x250	360	44
15.3963.01	20	180x320	450	46

U.N. APPROVED STACKABLE TANKS

Material: High density Polyethylene of high molecular weight
 Colour: Neutral - For food use
 They are equipped with self-sealing cap with total gasket in Polyethylene foam which ensure the tightness of liquids. Inner-cap for food. They can be supplied with stopcock PL 1394.

CONVENIENCE:

- The dimensions allow the overlap in the pallet 800x1200 and 1000x1200
- The handle allows for a secure and comfortable grip, the upper shape provides excellent stability
- The surface of the walls allows a safe and perfect application of labels
- The bottom has a concavity that allows an easy transfer

CERTIFICATIONS:

- They can be used for shipments of dangerous goods into the U.N. rules and in accordance with the regulation for road transport (ADR-RID), sea (IMO-IMDG), air (IATA-ICAO)

Item	Capacity lt	Dimension mm	Height mm	Mouth	U.N. Omologation
15.4110.01	5	163x187	236	45	3H1/Y1,8/150
15.4112.03	10	187x227	320	45	3H1/Y1,8/150
15.4114.02	20	236x288	380	45	3H1/Y1,8/180
15.4121.01	25	245x295	468	45	3H1/Y1,8/180
15.4118.01	30	287x314	446	45	3H1/Y1,8/180



WIDE MOUTH DRUMS

Material: High density Polyethylene - Color: Neutral - For food use
 Very strong and durable for extended use over time with carrying handle. With very large mouth to improve the phases of filling and emptying. On request we can provide stopcock BSP 3/4" threaded with a depth of 12 mm. with gasket for a perfect seal. They can be used with temperatures below 0°C. Cap with holes for warranty seal.



STOPCOCK 15.1393.02 TO BE ORDERED SEPARATELY

Material: High density Polyethylene 3/4" threaded
 For food use.
 Outlet hole ø 9 mm.

Item	Capacity lt	Diameter mm	Height mm	Mouth mm
15.1389.01	220	600	1030	350



WIDE MOUTH BOTTLE WITH HANDLE

Material: High density Polyethylene
 Color: Neutral - For food use

Item	Capacity ml	Dimension mm	Height mm	Mouth mm
15.3940.01	1000	85x105	220	45





HUBERLAB.

committed to science

Funnels - Vials Graduated jugs



Funnels



FUNNELS FOR TRANSFER

Material: High density Polyethylene

Color: White - For food use

General purpose for liquids. With rim, rib outer airvent and eyelet for easy hanging.

Item	Capacity ml	Upper diameter mm	Lower diameter mm	Tube lenght mm
15.1171.01	100	80	8	30
15.1172.01	300	100	10	40
15.1173.01	450	120	12	50
15.1174.01	750	140	16	55
15.1175.02	1500	180	18	67
15.1176.01	1750	210	21	80
15.1177.01	4000	260	25	85
15.1178.01	5500	310	30	90
15.1198.01	10000	420	40	125
15.1199.01	11000	450	45	150



FUNNEL WITH FILTER WITH FLEXIBLE SHANK

Material: High density Polyethylene

Color: White - For food use

Item	Upper diameter mm	Lower diameter mm
15.1202.09	210	20
15.1202.10	Spare filter	



FUNNELS FOR POWDERS

Material: Polypropylene

For food use

Designed for powder, they can be also used for liquids.

External ribs and high chemical resistance.

Item	Diameter mm	Tube diameter mm
15.1140.01	65	15
15.1141.01	80	21
15.1142.01	100	22
15.1143.01	120	26
15.1145.01	150	28

Funnels - Vials

FUNNELS FOR ANALYSIS

Material: Polypropylene

For food use - For standard laboratory uses

Autoclavable at + 120° C. With grooves to prevent air pockets and backpressure.

Item	Diameter mm	Tube lenght mm	Tube diameter mm
15.1212.01	40	40	6
15.1213.03	60	60	8
15.1214.01	80	80	8
15.1215.03	100	100	10
15.1216.01	120	120	16
15.1217.01	150	150	16



FUNNELS FOR ANALYSIS WITH INTERNAL RIBS TWIST

Material : Polypropylene

For food use - For standard laboratory uses.

Autoclavable at 120° C. Internal ribs that increase the filtration speed.

Item	Diameter mm	Tube lenght mm	Tube diameter mm
15.1212.02	40	40	6
15.1213.01	60	60	8
15.1214.02	80	80	8
15.1215.01	100	100	10
15.1216.02	120	120	16
15.1217.02	150	150	16



SCINTILLATION VIALS

Material: High density Polyethylene

Vials compatible with all liquid scintillation counters.

Leakproof closure with screw cap. heavy walled to minimize solvent losses.

Item	Capacity ml	Diameter mm	Diameter neck mm	Height mm	Cap size
15.9339.01	20	26,5	17,5	57,5	22,5



Graduated jugs



ROUND GRADUATED JUGS SHORT TYPE

Material: Polypropylene

For food use - Autoclavable

They resist to temperatures from -10°C to +120° C.

Color: Neutral - Very transparent especially when they contain liquid products.

Graduation is in embossed form, very strong with drip spout.

Item	Capacity ml	Graduation ml	Diameter mm	Height mm
15.1161.01	250	10	65	95
15.1156.01	500	25	75	120
15.1157.01	1000	50	120	135
15.1158.01	2000	100	145	170
15.1159.01	3000	100	165	185
15.1160.01	5000	250	190	235
• 15.1162.01	3000	500	165	185
• 15.1162.05	5000	500	190	235
15.1157.05	Lid with perforated spout			

• "SOFT TYPE" Material: Low density Polyethylene - Color: White - For food use



ROUND GRADUATED JUG TALL TYPE

Material: Polypropylene

For food use - Autoclavable

They resist to temperatures from -10 to +120° C

Color: Neutral - Very transparent especially when they contain liquid products. Graduation is in embossed form. Very strong with drip spout.

Item	Capacity ml	Diameter mm	Height mm
15.1092.02	250	75	115
15.1093.02	500	100	137
15.1094.04	1000	125	167
15.1095.02	2000	150	217
15.1096.04	3000	170	244



ROUND GRADUATED JUG TALL TYPE WITH BLUE GRADUATION

Material: Polypropylene

For food use - Autoclavable

They resist to temperatures from -10°C to +120°C.

Color: Neutral - Made with ultra transparent PP.

Thick walls and high chemical resistance. Handles with ergonomic shape to save space. Graduations is in blue enamel.

Item	Capacity ml	Diameter mm	Height mm
15.1092.03	250	75	115
15.1093.03	500	100	137
15.1094.05	1000	125	167
15.1095.03	2000	150	217
15.1096.05	3000	170	244

Cylinders

Beakers - Flasks



Graduated measuring cylinders



TALL TYPE CYLINDERS

Material: Polypropylene

For food use - Autoclavable at +121°C for 20 minutes

Made according to standards DIN 12681 and ISO 6706.

Graduation is in embossed form. They offer high chemical resistance; high stability base.

Item	Capacity ml	Graduation ml	Subdivision ml
15.1865.02	10	2	0,2
15.1077.03	25	5	0,5
15.1078.03	50	10	1,0
15.1079.03	100	10	1,0
15.1080.03	250	20	2,0
15.1081.04	500	50	5,0
15.1082.04	1000	100	10,0
15.1872.01	2000	200	20,0



TALL TYPE CYLINDERS WITH BLUE GRADUATION

Material: Polypropylene

For food use - Autoclavable at +121°C for 20 minutes

Made according to standards DIN 12681 and ISO 6706.

Graduation is blue enamel.

They offer high chemical resistance, high stability base.

Item	Capacity ml	Graduation ml	Subdivision ml
15.1865.01	10	2	0,2
15.1077.02	25	5	0,5
15.1078.02	50	10	1,0
15.1079.02	100	10	1,0
15.1080.02	250	20	2,0
15.1081.03	500	50	5,0
15.1082.03	1000	100	10,0



TALL TYPE CYLINDERS

Material: TPX

For food use - Transparent, autoclavable up to +170°C

Made according to standards ISO 6706.

Graduation in embossed form.

They offer high chemical resistance; high stability base.

Item	Capacity ml	Graduation ml	Subdivision ml
15.1340.01	10	2	0,2
15.1571.01	25	5	0,5
15.1572.01	50	10	1,0
15.1573.01	100	10	1,0
15.1574.01	250	20	2,0
15.1575.01	500	50	5,0
15.1576.01	1000	100	10,0
15.1347.01	2000	200	20,0

Graduated measuring cylinders - Beakers


A - SHORT TYPE CYLINDERS

Material: Polypropylene - For food use - Resistant up to 120°C

B - SHORT TYPE CYLINDERS

Material: TPX - For food use - Resistant up to 170°C

Item A	Capacity ml	Graduation ml	Subdivision ml	Item B
15.1814.01	10	2	0,25	15.1332.02
15.1815.01	25	5	1,0	15.1333.02
15.1816.01	50	10	2,5	15.1334.01
15.1817.01	100	25	5,0	15.1335.01
15.1820.01	250	50	10,0	15.1336.01
15.1821.01	500	100	10,0	15.1337.01
15.1822.01	1000	200	25,0	15.1338.01
15.1823.01	2000	400	50,0	15.1339.01



A - BEAKERS

Material: Polypropylene - For food use

Resistant up to 120°C

Autoclavable at +121°C for 20 minutes. Made according to standards ISO 7056. They offer high chemical resistance, high stability base.


B - BEAKERS WITH BLUE GRADUATION

Material: Polypropylene - For food use

Resistant up to 120°C

Graduation printed with blue enamel. Made according to standards ISO 7056.

Item A	Capacity ml	Subdivision ml	Item B
15.1801.02	25	1	15.1801.01
15.1802.02	50	2	15.1802.01
15.1803.02	100	5	15.1803.01
15.1805.02	250	10	15.1805.01
15.1806.02	500	10	15.1806.01
15.1808.02	1000	20	15.1808.01
15.1810.02	2000	50	15.1810.01
15.1811.01	3000	500	-
15.1812.01	5000	500	-



BEAKERS

Material: TPX

For food use

Transparent, autoclavable, heat resistance + 170°C

Item	Capacity ml	Subdivision ml
15.1321.24	25	1
15.1542.01	50	2
15.1543.01	100	5
15.1545.01	250	10
15.1546.01	500	10
15.1548.01	1000	20
15.1329.01	2000	50
15.1330.01	3000	500
15.1331.02	5000	500



Flasks



CONICAL FLASKS ERLENMEYER WITH SCREW CAP

Material: Polypropylene
For food and for laboratory use. Graduated.
They have standard neck and they are autoclavable.

Item	Capacity ml	Cap
15.1460.01	50	40
15.1461.01	100	40
15.1463.01	250	52
15.1464.02	500	52
15.1465.01	1000	52



VOLUMETRIC FLASKS WITH CAP

Material: Polypropylene
For food use
Transparent, high chemical resistance. Individually calibrated, they respect tolerances DIN/ISO class "B". Marked with blue enamel.
They are supplied with cap. They are autoclavable.

Item	Capacity ml	Height mm	Cap NS
15.1480.08	25	110	10/19
15.1481.07	50	140	12/21
15.1482.05	100	170	14/23
15.1483.01	250	220	14/23
15.1485.05	500	260	19/26
15.1486.01	1000	300	24/29



CONICAL CAPS


Material: High density Polyethylene
For food use
Conical caps as per DIN 12242. The body of the cap presents embossed rings, to ensure seal.

Item	Cap NS
15.4801.09	10/19
15.4811.01	12/21
15.4831.01	14/23
15.4851.01	19/26
15.4861.01	24/29

Tubing - Monobulb
Ballons - Stoppers
Stirring bars



Tubing



PVC TUBING


Material: PVC

Good flexibility and chemical resistance gives to this tubing a good use with apparatus assembled.

Autoclavable (121° C.)

Item	Inner Ø mm	Outer Ø mm	Pack mt
15.8110.01	3	5	20
15.8110.02	4	6	20
15.8110.03	5	8	20
15.8110.04	6	9	20

Item	Inner Ø mm	Outer Ø mm	Pack mt
15.8110.05	8	12	20
15.8110.06	10	14	20
15.8110.07	12	17	20



SILICON TUBING

Material: Silicon

Autoclavable (121° C.) - Translucent

Item	Inner Ø mm	Outer Ø mm	Pack mt
15.8120.01	1	3	10
15.8120.02	2	4	10
15.8120.03	3	5	10
15.8120.04	4	6	10

Item	Inner Ø mm	Outer Ø mm	Pack mt
15.8120.05	5	8	10
15.8120.06	6	9	10
15.8120.07	7	11	10
15.8120.08	8	12	10

Item	Inner Ø mm	Outer Ø mm	Pack mt
15.8120.09	9	13	10
15.8120.10	10	14	10
15.8121.03	12	17	10
15.8121.04	14	19	10



PTFE TUBING

Material: PTFE


Excellent resistance to acids, bases and solvents.

Temperature resistance from -200°C to +260°C.

Item	Inner Ø mm	Outer Ø mm	Pack mt
15.8130.01	2	4	5
15.8130.02	3	5	5
15.8130.03	4	6	5

Item	Inner Ø mm	Outer Ø mm	Pack mt
15.8130.04	6	8	5
15.8130.05	7	9	5
15.8130.06	8	10	5

Item	Inner Ø mm	Outer Ø mm	Pack mt
15.8130.07	9	11	5
15.8130.08	10	12	5
15.8130.09	12	14	5



RED RUBBER TUBING

Material: Red Rubber


Useful for transferring gas or non corrosive liquids.

Temperature resistance from -35°C to +70°C.

Item	Inner Ø mm	Outer Ø mm	Pack mt
15.8140.01	3	5	10
15.8140.02	4	7	10
15.8140.03	5	8	10

Item	Inner Ø mm	Outer Ø mm	Pack mt
15.8140.04	6	9	10
15.8140.05	6	12	10
15.8140.06	8	12	10

Item	Inner Ø mm	Outer Ø mm	Pack mt
15.8140.07	9	13	10
15.8140.08	10	14	10
15.8140.09	12	17	10



LATEX TUBING

Material: Latex

Used for transferring low density liquids.

Item	Inner Ø mm	Outer Ø mm	Pack mt
15.8150.01	4	7	10
15.8150.02	6	9	10
15.8150.03	8	12	10
15.8150.04	12	17	10


PARAVACUUM TUBING

Material: ParaRubber

Used for transferring air, acids, alkalines, solvents, basics.

Temperature resistance from -40°C to +70°C.

Item	Inner Ø mm	Outer Ø mm	Pack mt
15.8160.01	5	15	6
15.8160.02	6	18	6
15.8160.03	8	16	6
15.8160.04	12	21	6




Bulbs - Balloons - Stoppers

MONO BULB WITH TUBE

Material: Red Rubber
Pressure valves size 3


Item	Ø mm
15.8071.01	46



VACUUM BULB WITH TWO VALVES AND TUBE

Material: Red Rubber
Complete with vacuum and pressure valves.


Item	Ø mm
15.8072.01	46



RICHARDSON BULBS FOR AUTOMATIC BURETTES

Material: Red Rubber
With a net around the bulb.

Item	Ø mm
15.8073.09	46
15.8073.10	53




GAS SAMPLING BALLOONS

Material: Latex
Volume from 1000 to 8000 cc.

Item	Vol. cc	Ø mm
15.8060.09	1000	130
15.8060.10	1500	150
15.8060.11	2000	160

Item	Vol. cc	Ø mm
15.8060.13	4000	200
15.8060.14	5000	230
15.8060.15	8000	250



REVERSIBLE STOPPERS

Material: Silicon

Item	H/h mm	Ø inf. mm
15.8050.01	9/8	7
15.8050.02	14/8	10
15.8050.03	20/10	13
15.8050.04	20/10	15

Item	H/h mm	Ø inf. mm
15.8050.05	20/11	16
15.8050.06	18/12	20
15.8050.07	25/17	24
15.8050.08	28/22	31



Stoppers - Stirring bars - Latex bulb



SYLICON CONICAL STOPPERS

Material: Sylicon

Item	Ø sup. mm	Ø inf. mm	Height mm
15.8020.01	9	6	18
15.8020.02	12	9	18
15.8020.03	13	10	18
15.8020.04	14	11	18
15.8020.05	16	12	20
15.8020.06	18	13	22

Item	Ø sup. mm	Ø inf. mm	Height mm
15.8020.07	20	14	24
15.8020.08	23	16	26
15.8020.09	26	19	28
15.8020.10	30	21	29
15.8021.06	32	25	34
15.8021.07	37	28	37

Item	Ø sup. mm	Ø inf. mm	Height mm
15.8021.08	42	32	41
15.8021.09	45	34	45
15.8021.10	48	37	50
15.8021.11	54	41	53
15.8021.12	60	43	60



RED RUBBER CONICAL STOPPERS

Material: RedRubber

Item	Ø sup. mm	Ø inf. mm	Height mm
15.8030.01	9	6	18
15.8030.02	12	9	18
15.8030.03	13	10	18
15.8030.04	14	11	18
15.8030.05	16	12	20
15.8030.06	18	13	22

Item	Ø sup. mm	Ø inf. mm	Height mm
15.8030.07	20	14	24
15.8030.08	23	16	26
15.8030.09	26	19	28
15.8030.10	30	21	29
15.8031.01	32	25	34
15.8031.02	37	28	37

Item	Ø sup. mm	Ø inf. mm	Height mm
15.8031.03	42	32	41
15.8031.04	45	34	45
15.8031.05	48	37	50
15.8031.06	54	41	53
15.8031.07	60	43	60
15.8031.08	65	52	70




RED RUBBER CONICAL STOPPERS WITH 1 OR 2 HOLES

Material: RedRubber

Item 1 hole	Item 2 holes	Ø sup.	Ø inf.	Ø hole	Height mm
15.8040.03	-	12	9	3	18
15.8040.04	-	13	10	3	18
15.8040.05	-	14	11	3	18
15.8040.06	-	16	12	3	20
15.8040.07	15.8045.06	18	13	3	22
15.8040.08	15.8045.07	20	14	3	24
15.8040.09	15.8045.08	23	16	6	26

Item 1 hole	Item 2 holes	Ø sup. mm	Ø inf. mm	Ø hole mm	Height mm
15.8040.10	15.8045.09	26	19	6	28
15.8040.11	15.8045.10	30	21	6	29
15.8040.12	15.8045.11	32	25	6	34
15.8041.05	15.8045.12	37	28	6	37
15.8041.06	15.8045.13	42	32	6	41
15.8041.07	15.8045.14	45	34	6	45
15.8041.08	15.8045.15	48	37	6	50
15.8041.09	15.8046.10	54	41	6	53
15.8041.10	15.8046.11	60	43	6	60
15.8041.11	15.8046.12	65	52	6	70



MAGNETIC STIRRING BARS


Material: PTFE

Round smooth surface, for universal application.

Item	Ø mm	Length mm
15.8010.01	6	10
15.8010.02	6	15
15.8010.03	6	20
15.8010.04	6	25

Item	Ø mm	Length mm
15.8010.05	6	30
15.8010.06	6	35
15.8010.07	8	40
15.8010.08	8	45

Item	Ø mm	Length mm
15.8010.09	8	50
15.8010.10	9	60
15.8011.01	9	70



LATEX BULB FOR GLASS PASTEUR PIPETTES

Material: Latex

To be used with ranvier dropping tubes and pasteur pipettes.

Item	Height mm
15.8075.01	40

Containers Scoops



Stackable crates

HIGH DENSITY POLYETHYLENE CRATES WITH SEVERAL DIMENSIONS, SUITABLE TO CONTENT AND TRANSPORT EVERY KIND OF PRODUCTS, EASY TO CLEAN UNBREAKABLES - MATERIAL: HIGH DENSITY POLYETHYLENE - FOR FOOD USE CLOSED BOTTOM AND WALLS"



Item	Capacity lt	Weight gr	External measures mm	Internal measures mm	Color
15.3001.01	10	630	430x320x 95	390x290x 90	White
15.3003.01	16	1070	540x350x110	515x325x100	White
15.3005.02	16	1000	540x350x110	515x325x100	White

- with side holes for water drainage



Item	Capacity lt	Weight gr	External measures mm	Internal measures mm	Color
15.3004.02	55	2860	760x440x205	700x400x200	White



Item	Capacity lt	Weight gr	External measures mm	Internal measures mm	Color
15.1602.06	20	1070	440x315x200	380x280x195	White

Stackable crates

HIGH DENSITY POLYETHYLENE CRATES WITH SEVERAL DIMENSIONS, SUITABLE TO CONTENT AND TRANSPORT EVERY KIND OF PRODUCTS, EASY TO CLEAN
UNBREAKABLES - MATERIAL: HIGH DENSITY POLYETHYLENE - FOR FOOD USE CLOSED BOTTOM AND WALLS"

Item	Capacity lt	Weight gr	External measures mm	Internal measures mm	Color
15.1603.01	50	2150	625x365x285	560x330x280	White



Item	Capacity lt	Weight gr	External measures mm	Internal measures mm	Color
15.1604.06	80	3170	725x440x330	650x380x320	White



Item	Capacity lt	Weight gr	External measures mm	Internal measures mm	Color
15.1601.09	100	4300	600x500x445	535x450x430	White



Ice box - Scoop - Large cup




ISOTHERMAL CONTAINERS

Material: High density Polyethylene

Color: Blue - For food use

The thermal insulation is ensured by the injection of Polyurethane between the walls of the container. Suitable for transporting substances where it is necessary to maintain a constant temperature.

Item	Capacity lt	Dimension mm	Height mm
15.1780.01	20	225x420	330
15.1781.01	25	250x420	350
15.1782.01	32	270x440	380
15.1783.01	45	305x520	400



A - SCOOP WITH FLAT BOTTOM

Material: Polypropylene

Color: Neutral - For food use

For withdrawal of solids and granulate samples.

The measure includes with handle.

B - LARGE CUP

Material: High density Polyethylene

Color: White - For food use

Item	Capacity ml	Total length mm	Useful length mm
15.1066.01	2000	400	305

Item	Capacity ml	Total length mm	Useful length mm
15.2147.01	3000	200	115



SCOOPS

Material: High density Polyethylene

Color: Neutral - For food use

For withdrawal of solid and granulates samples.

They have the rounded bottom for a better use.

Item	Capacity ml	Total length mm	Useful length mm
15.1062.02	90	180	125
15.1063.02	250	240	155
15.1064.03	600	300	195
15.1065.01	1100	360	240
15.1067.02	2500	400	285

Spatulas - Pump

SPATULAS

Material: High density Polyethylene

Color: White - For food use

Very strong. Suitable for the mixing of chemicals and food.

Item	Length mm
15.2032.01	600
15.2033.01	800
15.2035.01	1000
15.2038.01	1200
15.2039.01	1400



PUMP FOR TRANSFER

Material: Polypropylene

Piston operation, telescopic suction pipe. Solid structure, suitable for heavy-duty applications. 2" threaded port.

Item	Capacity ml
15.1955.01	250




PUMP

Material: High density Polyethylene/PVC

For food use

Item	
15.1935.01	



Bottles carrier - Bed basin - Maked bed bottle



BOTTLES CARRIER

Material: High density Polyethylene

Color: Red

Designed for 6 bottles with a diameter up to 95 mm.

Produced in hdpe with an ergonomic design of handles.

Item	Dimension mm
15.1402.03	320x300x200



BED BASIN

Material: Polypropylene

Color: Neutral

Produced in polypropylene easy to clean, it can resist up to 120° C.

Item	
15.1401.06	



MALE BED BOTTLE

Material: High density Polyethylene

Color: Neutral

Item	
15.1400.08	

- independent, family-owned Swiss Company
- specialized retailer
- established in 1968
- comprehensive productportfolio for laboratories consumables and instruments
- largest warehouse in Switzerland for laboratory consumables
- own delivery service
- eShop, e-commerce & B2B



Price List

Page	Item	Description	Pieces per unit	Price per Unit CHF
2	15.1408.01	Wide mouth bottle 50 ml, graduated, LDPE	1	0.60
2	15.1409.01	Wide mouth bottle 100 ml, graduated, LDPE	1	0.70
2	15.1410.03	Wide mouth bottle 250 ml, graduated, LDPE	1	1.00
2	15.1411.01	Wide mouth bottle 500 ml, graduated, LDPE	1	1.05
2	15.1412.01	Wide mouth bottle 1000 ml, graduated, LDPE	1	1.85
2	15.1413.01	Wide mouth bottle 2000 ml, graduated, LDPE	1	3.45
2	15.1608.06	Square bottle, wide mouth 25 ml, graduated, HDPE	1	0.60
2	15.1609.01	Square bottle, wide mouth 50 ml, graduated, HDPE	1	0.70
2	15.1610.08	Square bottle, wide mouth 100 ml, graduated, HDPE	1	0.70
2	15.1611.01	Square bottle, wide mouth 250 ml, graduated, HDPE	1	0.95
2	15.1612.06	Square bottle, wide mouth 500 ml, graduated, HDPE	1	1.30
2	15.1613.01	Square bottle, wide mouth 1000 ml, graduated, HDPE	1	2.00
2	15.1614.02	Square bottle, wide mouth 2000 ml, graduated, HDPE	1	3.50
3	15.1317.01	Narrow mouth bottle 50 ml, graduated, LDPE	1	0.65
3	15.1319.01	Narrow mouth bottle 125 ml, graduated, LDPE	1	0.75
3	15.1323.02	Narrow mouth bottle 250 ml, graduated, LDPE	1	0.95
3	15.1324.02	Narrow mouth bottle 500 ml, graduated, LDPE	1	1.15
3	15.1325.02	Narrow mouth bottle 1000 ml, graduated, LDPE	1	1.85
3	15.1326.02	Narrow mouth bottle 2000 ml, graduated, LDPE	1	3.65
3	15.1323.03	Narrow mouth bottle 250 ml, LDPE, grey	1	1.15
3	15.1324.03	Narrow mouth bottle 500 ml, LDPE, grey	1	1.35
3	15.1325.03	Narrow mouth bottle, 1000 ml, LDPE, grey	1	2.30
3	15.1326.03	Narrow mouth bottle, 2000 ml, LDPE, grey	1	4.40
4	15.9126.01	GL 32 Wide mouth bottle PP, 50 ml - with cap	1	0.95
4	15.9326.01	GL 32 Wide mouth bottle PP, 100 ml - with cap	1	1.25
4	15.9526.01	GL 40 Wide mouth bottle PP, 250 ml - with cap	1	1.60
4	15.9726.01	GL 50 Wide mouth bottle PP, 500 ml - with cap	1	2.05
4	15.9926.01	GL 65 Wide mouth bottle PP, 1000 ml - with cap	1	4.30
4	15.9136.01	GL 65 Wide mouth bottle PP, 2000 ml - with cap	1	6.65
4	15.9806.01	GL 32 Wide mouth bottle PE 50 ml with cap	1	0.85
4	15.9016.01	GL 32 Wide mouth bottle PE 100 ml with cap	1	1.10
4	15.9216.01	GL 40 Wide mouth bottle PE 250 ml with cap	1	1.45
4	15.9416.01	GL 50 Wide mouth bottle PE 500 ml with cap	1	1.85
4	15.9616.01	GL 65 Wide mouth bottle PE 1000 ml with cap	1	3.95
4	15.9816.01	GL 65 Wide mouth bottle PE 2000 ml with cap	1	6.10
5	15.3608.01	GL 32 Wide mouth bottle PP 50 ml with cap, graduated	1	1.40
5	15.3610.01	GL 32 Wide mouth bottle PP 100 ml with cap, graduated	1	1.70
5	15.3612.01	GL 40 Wide mouth bottle PP 250 ml with cap, graduated	1	2.10
5	15.3614.01	GL 50 Wide mouth bottle PP 500 ml with cap, graduated	1	2.30
5	15.3616.01	GL 65 Wide mouth bottle PP 1000 ml with cap, graduated	1	4.80
5	15.3618.01	GL 65 Wide mouth bottle PP 2000 ml with cap, graduated	1	8.50
5	15.2608.01	GL 32 Wide mouth bottle PE, amber 50 ml with cap	1	1.35
5	15.2610.01	GL 32 Wide mouth bottle PE, amber 100 ml with cap	1	1.90
5	15.2612.01	GL 40 Wide mouth bottle PE, amber 250 ml with cap	1	2.30
5	15.2614.01	GL 50 Wide mouth bottle PE, amber 500 ml with cap	1	2.90
5	15.2616.01	GL 65 Wide mouth bottle PE, amber 1000 ml with cap	1	6.20
6	15.9595.01	GL 14 Narrow mouth bottle PP, 20 ml - with cap	1	0.50
6	15.9795.01	GL 18 Narrow mouth bottle PP, 100 ml - with cap	1	0.90
6	15.9995.01	GL 25 Narrow mouth bottle PP, 250 ml - with cap	1	1.40
6	15.9106.01	GL 25 Narrow mouth bottle PP, 500 ml - with cap	1	1.70
6	15.9306.01	GL 28 Narrow mouth bottle PP, 1000 ml - with cap	1	3.15
6	15.2620.01	GL 18 Narrow mouth bottle PE 100 ml with cap	1	1.30
6	15.2622.01	GL 25 Narrow mouth bottle PE 250 ml with cap	1	2.00
6	15.2624.01	GL 25 Narrow mouth bottle PE 500 ml with cap	1	2.40
6	15.2626.01	GL 28 Narrow mouth bottle PE 1000 ml with cap	1	4.50
7	15.9108.01	GL 18 Narrow mouth bottle PE, 50 ml without cap	1	0.55
7	15.9167.95	GL 18 Narrow mouth bottle PE, 100 ml without cap	1	0.65
7	15.9167.99	GL 25 Narrow mouth bottle PE, 250 ml without cap	1	1.05
7	15.9168.01	GL 25 Narrow mouth bottle PE, 500 ml without cap	1	1.35
7	15.9188.01	GL 28 Narrow mouth bottle PE, 1000 ml without cap	1	2.60
7	15.9109.02	GL 28 Narrow mouth bottle PE, 2000 ml without cap	1	4.80
7	15.9008.01	GL 18 cap HDPE	1	0.25
7	15.9048.01	GL 25 cap HDPE	1	0.35
7	15.9088.01	GL 28 cap HDPE	1	0.35
7	15.9603.01	GL 18 long dropper	1	0.40
7	15.9703.01	GL 25 long dropper	1	0.45
7	15.9803.01	GL 28 long dropper	1	0.50
7	15.9303.01	GL 18 short dropper	1	0.40
7	15.9403.01	GL 25 short dropper	1	0.45
7	15.9003.01	GL 18 screw cap for PE wash bottle	1	1.25
7	15.9103.01	GL 25 screw cap for PE wash bottle	1	1.55
7	15.9203.01	GL 28 screw cap for PE wash bottle	1	1.70
8	15.3710.01	Wide neck bottle, PP 100 ml without exagonal cap 29/32	1	4.00

Page	Item	Description	Pieces per unit	Price per Unit CHF
8	15.3714.01	Wide neck bottle, PP 500 ml without exagonal cap 45/40	1	6.10
8	15.3716.01	Wide neck bottle, PP 1000 ml without exagonal cap 60/46	1	8.25
8	15.3720.01	Narrow neck bottle, PP 100 ml without exagonal cap 14/23	1	2.70
8	15.3722.01	Narrow neck bottle, PP 250 ml without exagonal cap 19/26	1	3.45
8	15.3724.01	Narrow neck bottle, PP 500 ml without exagonal cap 24/29	1	4.65
8	15.3726.01	Narrow neck bottle, PP 1000 ml without exagonal cap 29/32	1	6.00
8	15.2240.07	Hexagonal conical cap PE, 7/16	50	20.00
8	15.2241.03	Hexagonal conical cap PE, 10/19	50	27.50
8	15.2242.01	Hexagonal conical cap PE, 12/21	50	32.50
8	15.2243.01	Hexagonal conical cap PE, 14/23	50	35.00
8	15.2244.01	Hexagonal conical cap PE, 14/15	50	35.00
8	15.2245.01	Hexagonal conical cap PE, 18/26	50	40.00
8	15.2246.01	Hexagonal conical cap PE, 16/26	50	40.00
8	15.2247.01	Hexagonal conical cap PE, 24/29	50	47.50
8	15.2248.01	Hexagonal conical cap PE, 24/20	50	47.50
8	15.2249.01	Hexagonal conical cap PE, 29/32	50	60.00
8	15.2250.01	Hexagonal conical cap PE, 29/22	50	60.00
8	15.2251.01	Hexagonal conical cap PE, 34/35	50	112.50
8	15.2252.01	Hexagonal conical cap PE, 45/40	50	134.50
8	15.2253.01	Hexagonal conical cap PE, 60/46	50	182.50
8	15.2254.01	Hexagonal conical cap PE, 85/65	50	255.00
9	15.1435.01	Round bottle, HDPE, with original seal cap, Ø 41 mm (DIN 50), 1000 ml	1	1.90
9	15.1451.02	Approved U.N. bottles, HDPE, 250 ml	1	1.95
9	15.1452.01	Approved U.N. bottles, HDPE, 500 ml	1	2.15
9	15.1453.02	Approved U.N. bottles, HDPE, 1000 ml	1	2.60
9	15.1351.02	HDPE narrow mouth bottle with original seal cap, 50 ml	1	0.60
9	15.1352.02	HDPE narrow mouth bottle with original seal cap, 100 ml	1	0.60
9	15.1352.03	HDPE narrow mouth bottle with original seal cap, 100 ml (16'000 pcs)	1	0.65
9	15.1353.02	HDPE narrow mouth bottle with original seal cap, 250 ml	1	0.85
9	15.1354.03	HDPE narrow mouth bottle with original seal cap, 500 ml	1	1.05
9	15.1355.01	HDPE narrow mouth bottle with original seal cap, 750 ml	1	1.35
9	15.1356.01	HDPE narrow mouth bottle with original seal cap, 1000 ml	1	1.50
10	15.1343.02	Wide mouth bottles Ø 30 mm, with original seal cap, 250 ml	1	1.10
10	15.1344.01	Wide mouth bottles Ø 30 mm, with original seal cap, 500 ml	1	1.35
10	15.1345.02	Wide mouth bottles Ø 30 mm, with original seal cap, 750 ml	1	1.60
10	15.1346.01	Wide mouth bottles Ø 30 mm, with original seal cap, 1000 ml	1	1.80
10	15.1373.01	HDPE bottle with child proof cap, 250 ml	1	1.10
10	15.1374.01	HDPE bottle with child proof cap, 500 ml	1	0.95
10	15.1376.01	HDPE bottle with child proof cap, 1000 ml	1	1.45
10	15.1363.01	HDPE drop bottle with nozzle, 250 ml	1	0.90
10	15.1364.01	HDPE drop bottle with nozzle, 500 ml	1	1.05
10	15.1366.01	HDPE drop bottle with nozzle, 1000 ml	1	1.50
11	15.1523.48	Square bottle HDPE, narrow mouth, graduated with original seal cap, 250 ml	1	1.15
11	15.1525.01	Square bottle HDPE, narrow mouth, graduated with original seal cap, 500 ml	1	1.45
11	15.1526.01	Square bottle HDPE, narrow mouth, graduated with original seal cap, 1000 ml	1	1.90
11	15.1524.01	Square bottle HDPE, wide mouth, graduated with original seal cap, 250 ml	1	1.20
11	15.1535.01	Square bottle HDPE, wide mouth, graduated with original seal cap, 500 ml	1	1.55
11	15.1536.01	Square bottle HDPE, wide mouth, graduated with original seal cap, 1000 ml	1	2.10
11	15.1380.02	HDPE narrow mouth bottle with original seal cap and inner-cap, 25 ml	1	0.70
11	15.1381.02	HDPE narrow mouth bottle with original seal cap and inner-cap, 50 ml	1	0.80
11	15.1382.02	HDPE narrow mouth bottle with original seal cap and inner-cap, 100 ml	1	0.85
11	15.1383.02	HDPE narrow mouth bottle with original seal cap and inner-cap, 250 ml	1	0.90
11	15.1384.04	HDPE narrow mouth bottle with original seal cap and inner-cap, 500 ml	1	1.05
11	15.1385.02	HDPE narrow mouth bottle with original seal cap and inner-cap, 1000 ml	1	1.50
12	15.1302.03	HDPE wide mouth bottles, short type, 100 ml	1	0.80
12	15.1303.03	HDPE wide mouth bottles, short type, 150 ml	1	0.85
12	15.1304.03	HDPE wide mouth bottles, short type, 200 ml	1	0.90

valid until 31.12.2017

Page	Item	Description	Pieces per unit	Price per Unit CHF
12	15.1036.01	HDPE wide mouth bottles with screw cap, 2000 ml	1	3.45
12	15.1037.03	HDPE wide mouth bottles with screw cap, 3000 ml	1	3.75
12	15.1394.03	Spray bottles LDPE, 500 ml	1	4.95
12	15.1395.02	Spray bottles LDPE, 900 ml	1	4.95
12	15.1396.03	Spray nozzle	1	2.80
13	15.1101.02	LDPE drop bottle, 50 ml	1	1.15
13	15.1102.03	LDPE drop bottle, 125 ml	1	2.25
13	15.1103.03	LDPE drop bottle, 250 ml	1	2.75
13	15.1105.01	LDPE drop bottle, 500 ml	1	3.30
13	15.1106.03	LDPE drop bottle, 750 ml	1	3.75
13	15.1107.06	LDPE drop bottle, 1000 ml	1	4.30
13	15.1180.01	Wash bottle LDPE, narrow neck 100 ml	1	2.25
13	15.1182.01	Wash bottle LDPE, narrow neck 250 ml	1	2.45
13	15.1183.01	Wash bottle LDPE, narrow neck 500 ml	1	2.90
13	15.1184.01	Wash bottle LDPE, narrow neck 750 ml	1	3.40
13	15.1185.01	Wash bottle LDPE, narrow neck 1000 ml	1	3.90
13	15.8900.01	LABEL DISTILLED WATER	1	2.25
13	15.8900.02	LABEL ACETONE	1	2.25
13	15.8900.03	LABEL METHANOL	1	2.25
13	15.8900.04	LABEL ISOPROPANOL	1	2.25
13	15.8900.05	LABEL ETHYL ACETATE	1	2.25
13	15.8900.06	LABEL ETHANOL	1	2.25
13	15.8900.07	LABEL ETHANOL	1	2.25
13	15.8900.08	LABEL N-HEXANE	1	2.25
13	15.8900.09	LABEL METH. CHLORIDE	1	2.25
13	15.8900.10	LABEL TOLUOL	1	2.25
13	15.8901.01	LABEL XYLOL	1	2.25
14	15.1470.08	Wash bottle wide neck ACETONE 500 ml	1	7.65
14	15.1471.01	Wash bottle wide neck DISTILLED WATER 500 ml	1	7.65
14	15.1472.08	Wash bottle wide neck ETHANOL	1	7.65
14	15.1473.01	Wash bottle wide neck ETHYL ACETATE	1	7.65
14	15.1474.01	Wash bottle wide neck ISOPROPANOL	1	7.65
14	15.1475.01	Wash bottle wide neck METHANOL	1	7.65
14	15.1476.01	Wash bottle wide neck METHYLETHYLKETONE	1	7.65
14	15.1477.01	Wash bottle wide neck TETRAHYDROFURAN	1	7.65
14	15.1478.01	Wash bottle wide neck TOLUENE	1	7.65
14	15.1031.01	HDPE Wide mouth bottles with screw cap, 2000 ml	1	12.85
14	15.1032.03	HDPE Wide mouth bottles with screw cap, 3000 ml	1	13.75
14	15.1448.01	HDPE Wide mouth bottles with original seal cap, 2000 ml	1	5.40
16	15.1020.05	HDPE jars with screw cap and inner-cap, 70 ml	1	0.80
16	15.1021.01	HDPE jars with screw cap and inner-cap, 120 ml	1	0.85
16	15.1022.01	HDPE jars with screw cap and inner-cap, 250 ml	1	0.90
16	15.1023.02	HDPE jars with screw cap and inner-cap, 500 ml	1	1.25
16	15.1025.10	HDPE jars with screw cap and inner-cap, 750 ml	1	1.50
16	15.1024.01	HDPE jars with screw cap and inner-cap, 1000 ml	1	1.85
16	15.1026.02	HDPE jars with screw cap and inner-cap, 2000 ml	1	2.40
16	15.1054.01	HDPE jars with original seal cap, 500 ml	1	1.60
16	15.1027.01	HDPE jars with original seal cap, 1000 ml	1	1.90
16	15.1028.02	HDPE jars with original seal cap, 1500 ml	1	2.25
16	15.1029.01	HDPE jars with original seal cap, 2000 ml	1	2.50
17	15.1043.02	HDPE square base jar, with pressure sealing cap, white, 700 ml	1	3.25
17	15.1045.02	HDPE square base jar, with pressure sealing cap, white, 1300 ml	1	3.85
17	15.1055.43	Stackable square base jar, with sealing cap, 5 ml, HDPE	1	15.40
17	15.1058.02	Stackable square base jar, with sealing cap, 10 ml, HDPE	1	19.65
17	15.1539.01	HDPE square base jar, with screw self-sealing cap, 1000 ml	1	3.70
18	15.2195.01	PP jar with original pressure seal cap, 650 ml	1	2.85
18	15.2415.01	Body of 50 ml jar, PP, white	1	0.40
18	15.2415.02	Pressure seal cap for 50 ml jar, white	1	0.15
18	15.2418.01	Body of 100 ml jar, PP, white	1	0.55
18	15.2418.02	Pressure seal cap for 100 ml jar, white	1	0.25
18	15.2420.01	Body of 150 ml jar, PP, white	1	0.80
18	15.2420.02	Pressure seal cap for 150 ml jar, white	1	0.30
18	15.2426.01	Body of 280 ml jar, PP, white	1	1.10
18	15.2426.02	Pressure seal cap for 280 ml jar, white	1	0.40
18	15.2430.01	Body of 500 ml jar, PP, white	1	1.65
18	15.2430.02	Pressure seal cap for 500 ml jar, white	1	0.75
18	15.2433.01	Body of 1000 ml jar, PP, white	1	2.75
18	15.2433.02	Pressure seal cap for 1000 ml jar, white	1	1.05
18	15.2434.01	Body of 1500 ml jar, PP, white	1	4.05
18	15.2434.02	Pressure seal cap for 1500 ml jar, white	1	1.05
18	15.2450.01	Body of 50 ml jar, PP, white	1	0.30
18	15.2450.02	Seal cap for 50 ml jar, white	1	0.30
18	15.2452.01	Body of 100 ml jar, PP, white	1	0.45
18	15.2452.02	Seal cap for 100 ml jar, white	1	0.30
18	15.2454.01	Body of 250 ml jar, PP, white	1	0.95

Page	Item	Description	Pieces per unit	Price per Unit CHF
18	15.2458.01	Body of 500 ml jar, PP, white	1	1.30
18	15.2458.02	Seal cap for 500 ml jar, white	1	0.85
18	15.2462.01	Body of 1000 ml jar, PP, white	1	2.40
18	15.2462.02	Seal cap for 1000 ml jar, white	1	1.25
18	15.2464.01	Body of 1250 ml jar, PP, white	1	2.95
18	15.2464.02	Seal cap for 1250 ml jar, white	1	1.25
19	15.1243.01	PP conical jar, 300 ml, with cap	1	1.10
19	15.1245.01	PP conical jar, 650 ml, with cap	1	1.40
19	15.1246.01	PP conical jar, 1000 ml, with cap	1	1.85
19	15.2174.01	30 ml transparent jar, PS	1	0.35
19	15.2175.05	50 ml transparent jar, PS	1	0.45
19	15.2177.01	100 ml transparent jar, PS	1	0.65
19	15.2178.01	150 ml transparent jar, PS	1	0.80
19	15.2180.01	200 ml transparent jar, PS	1	0.95
19	15.2182.01	250 ml transparent jar, PS	1	1.05
19	15.2184.01	500 ml transparent jar, PS	1	2.55
19	15.1307.01	HDPE graduated bottle with brush and steel ball, 50 ml	1	2.20
22	15.1187.01	HDPE carboy without stopcock, cap with 1 hole, 5 Lt	1	11.50
22	15.1188.01	HDPE carboy without stopcock, cap with 1 hole, 10 Lt	1	14.40
22	15.1190.50	HDPE carboy without stopcock, cap with 1 hole, 15 Lt	1	20.15
22	15.1190.55	HDPE carboy without stopcock, cap with 1 hole, 20 Lt	1	26.95
22	15.1191.01	HDPE carboy without stopcock, cap with 1 hole, 25 Lt	1	40.20
22	15.1192.01	HDPE carboy without stopcock, cap with 1 hole, 50 Lt	1	60.35
22	15.1187.02	HDPE carboy with stopcock, cap with 1 hole, 5 Lt	1	14.40
22	15.1188.02	HDPE carboy with stopcock, cap with 1 hole, 10 Lt	1	17.20
22	15.1190.00	HDPE carboy with stopcock, cap with 1 hole, 15 Lt	1	23.20
22	15.1190.02	HDPE carboy with stopcock, cap with 1 hole, 20 Lt	1	30.00
22	15.1191.02	HDPE carboy with stopcock, cap with 1 hole, 25 Lt	1	42.20
22	15.1192.02	HDPE carboy with stopcock, cap with 1 hole, 50 Lt	1	61.80
22	15.1393.02	3/4" stopcock	1	3.35
22	15.1386.02	HDPE wide mouth drum, cap with 1 hole, 30 Lt	1	47.30
22	15.1387.02	HDPE wide mouth drum, cap with 1 hole, 50 Lt	1	71.75
22	15.1388.02	HDPE wide mouth drum, cap with 1 hole, 100 Lt	1	110.20
22	15.1393.02	3/4" stopcock	1	3.35
23	15.1505.01	HDPE bottle with cap and inner cap, 5 lt	1	7.10
23	15.1510.06	HDPE bottle with cap and inner cap, 10 lt	1	8.30
23	15.1414.01	HDPE wide mouth bottle, cap with hole and inner cap, 5 lt	1	9.55
23	15.1415.01	HDPE wide mouth bottle, cap with hole and inner cap, 10 lt	1	13.00
23	15.1416.01	HDPE wide mouth bottle, cap with hole and inner cap, 15 lt	1	20.55
23	15.1417.01	HDPE wide mouth bottle, cap with hole and inner cap, 20 lt	1	23.80
23	15.1350.02	HDPE wide mouth bottle, cap with hole and inner cap, 25 lt	1	28.90
24	15.3941.01	HDPE jerrycans with original seal cap 1 lt	1	2.70
24	15.3942.01	HDPE jerrycans with original seal cap 2 lt	1	3.50
24	15.3943.01	HDPE jerrycans with cap 3 lt	1	5.10
24	15.3948.01	HDPE jerrycans with original seal cap 5 lt	1	5.15
24	15.3949.01	HDPE jerrycans with original seal cap 10 lt	1	7.45
24	15.3958.01	HDPE jerrycans with original seal cap 15 lt	1	10.75
24	15.3950.04	HDPE jerrycans with original seal cap 20 lt	1	11.85
24	15.3951.01	HDPE jerrycans with original seal cap 25 lt	1	15.85
24	15.3952.01	HDPE jerrycans with original seal cap 30 lt	1	17.45
24	15.3961.01	HDPE Jerrycans with screwed stopcock, 10 lt	1	11.95
24	15.3963.01	HDPE Jerrycans with screwed stopcock, 20 lt	1	25.25
24	15.5953.01	Plastirob Stopcock	1	2.05
25	15.4110.01	HDPE jerrycans with original seal cap, U.N. approved, with inner cap, 5 lt	1	4.95
25	15.4112.03	HDPE jerrycans with original seal cap, U.N. approved, with inner cap, 10 lt	1	7.55
25	15.4114.02	HDPE jerrycans with original seal cap, U.N. approved, with inner cap, 20 lt	1	14.50
25	15.4121.01	HDPE jerrycans with original seal cap, U.N. approved, with inner cap, 25 lt	1	15.85
25	15.4118.01	HDPE jerrycans with original seal cap, U.N. approved, with inner cap, 30 lt	1	18.00
25	15.1389.01	Drum of HDPE, 220 lt	1	214.70
25	15.1393.02	3/4" stopcock	1	3.35
25	15.3940.01	HDPE wide mouth bottle with handle, 1000 ml	1	3.35
28	15.1171.01	Funnel 8 cm Ø, HDPE	1	0.65
28	15.1172.01	Funnel 10 cm Ø, HDPE	1	1.30
28	15.1173.01	Funnel 12 cm Ø, HDPE	1	1.50
28	15.1174.01	Funnel 14 cm Ø, HDPE	1	1.95
28	15.1175.02	Funnel 17 cm Ø, HDPE	1	3.70
28	15.1176.01	Funnel 21 cm Ø, HDPE	1	4.80
28	15.1177.01	Funnel 25 cm Ø, HDPE	1	7.15
28	15.1178.01	Funnel 30 cm Ø, HDPE	1	9.55
28	15.1198.01	Funnel 40 cm Ø, HDPE	1	17.40
28	15.1199.01	Funnel 45 cm Ø, HDPE	1	19.65

Price List

Page	Item	Description	Pieces per unit	Price per Unit CHF
28	15.1140.01	Funnel for powder, PP, 65 mm	1	5.65
28	15.1141.01	Funnel for powder, PP, 80 mm	1	8.10
28	15.1142.01	Funnel for powder, PP, 100 mm	1	8.35
28	15.1143.01	Funnel for powder, PP, 120 mm	1	10.45
28	15.1145.01	Funnel for powder, PP, 150 mm	1	15.50
29	15.1212.01	PP funnel for analysis, Ø 40 mm	1	0.95
29	15.1213.03	PP funnel for analysis, Ø 60 mm	1	1.20
29	15.1214.01	PP funnel for analysis, Ø 80 mm	1	1.70
29	15.1215.03	PP funnel for analysis, Ø 100 mm	1	2.50
29	15.1216.01	PP funnel for analysis, Ø 120 mm	1	3.65
29	15.1217.01	PP funnel for analysis, Ø 150 mm	1	4.60
29	15.1212.02	PP funnel for analysis, Ø 40 mm, with rib	1	1.10
29	15.1214.02	PP funnel for analysis, Ø 80 mm, with rib	1	1.85
29	15.1215.01	PP funnel for analysis, Ø 100 mm, with rib	1	2.95
29	15.1216.02	PP funnel for analysis, Ø 120 mm, with rib	1	4.05
29	15.1217.02	PP funnel for analysis, Ø 150 mm, with rib	1	5.20
29	15.9339.01	Scintillation vials, HDPE, with screw cap, 20 ml	1000	250.00
30	15.1161.01	PP jug short type with handle and spout, graduated, 250 ml	1	2.50
30	15.1156.01	PP jug short type with handle and spout, graduated, 500 ml	1	1.75
30	15.1157.01	PP jug short type with handle and spout, graduated, 1000 ml	1	1.90
30	15.1158.01	PP jug short type with handle and spout, graduated, 2000 ml	1	3.70
30	15.1159.01	PP jug short type with handle and spout, graduated, 3000 ml	1	4.55
30	15.1160.01	PP jug short type with handle and spout, graduated, 5000 ml	1	7.05
30	15.1162.01	LDPE jug short type with handle and spout, graduated, white, 3000 ml	1	6.05
30	15.1162.05	LDPE jug short type with handle and spout, graduated, white, 5000 ml	1	8.15
30	15.1157.05	Lid for 1000 ml jug with perforated spout	1	0.65
30	15.1092.02	PP jug with handle and spout, graduated, 250 ml	1	4.55
30	15.1093.02	PP jug with handle and spout, graduated, 500 ml	1	6.30
30	15.1094.04	PP jug with handle and spout, graduated, 1000 ml	1	8.45
30	15.1095.02	PP jug with handle and spout, graduated, 2000 ml	1	14.35
30	15.1096.04	PP jug with handle and spout, graduated, 3000 ml	1	20.45
30	15.1092.03	PP jug with handle and spout, blue graduated, 250 ml	1	3.80
30	15.1093.03	PP jug with handle and spout, blue graduated, 500 ml	1	4.10
30	15.1094.05	PP jug with handle and spout, blue graduated, 1000 ml	1	6.25
30	15.1095.03	PP jug with handle and spout, blue graduated, 2000 ml	1	11.65
30	15.1096.05	PP jug with handle and spout, blue graduated, 3000 ml	1	16.25
32	15.1865.02	Tall type cylinders, PP, graduated, 10 ml	1	2.35
32	15.1077.03	Tall type cylinders, PP, graduated, 25 ml	1	2.95
32	15.1078.03	Tall type cylinders, PP, graduated, 50 ml	1	3.30
32	15.1079.03	Tall type cylinders, PP, graduated, 100 ml	1	4.30
32	15.1080.03	Tall type cylinders, PP, graduated, 250 ml	1	6.75
32	15.1081.04	Tall type cylinders, PP, graduated, 500 ml	1	10.25
32	15.1082.04	Tall type cylinders, PP, graduated, 1000 ml	1	16.70
32	15.1872.01	Tall type cylinders, PP, graduated, 2000 ml	1	47.65
32	15.1865.01	Tall type cylinders, PP, blue graduation, 10 ml	1	3.20
32	15.1077.02	Tall type cylinders, PP, blue graduation, 25 ml	1	3.75
32	15.1078.02	Tall type cylinders, PP, blue graduation, 50 ml	1	4.30
32	15.1079.02	Tall type cylinders, PP, blue graduation, 100 ml	1	5.65
32	15.1080.02	Tall type cylinders, PP, blue graduation, 250 ml	1	8.80
32	15.1081.03	Tall type cylinders, PP, blue graduation, 500 ml	1	12.95
32	15.1082.03	Tall type cylinders, PP, blue graduation, 1000 ml	1	21.85
32	15.1340.01	Tall type cylinders, TPX, graduated, 10 ml	1	6.95
32	15.1571.01	Tall type cylinders, TPX, graduated, 25 ml	1	9.55
32	15.1572.01	Tall type cylinders, TPX, graduated, 50 m	1	10.80
32	15.1573.01	Tall type cylinders, TPX, graduated, 100 ml	1	11.50
32	15.1574.01	Tall type cylinders, TPX, graduated, 250 ml	1	17.90
32	15.1575.01	Tall type cylinders, TPX, graduated, 500 ml	1	26.60
32	15.1576.01	Tall type cylinders, TPX, graduated, 1000 ml	1	41.65
32	15.1347.01	Tall type cylinders, TPX, graduated, 2000 ml	1	82.80
33	15.1814.01	Short type cylinders, PP, graduated, 10 ml	1	5.40
33	15.1815.01	Short type cylinders, PP, graduated, 25 ml	1	6.00
33	15.1816.01	Short type cylinders, PP, graduated, 50 ml	1	2.00
33	15.1817.01	Short type cylinders, PP, graduated, 100 ml	1	2.70
33	15.1820.01	Short type cylinders, PP, graduated, 250 ml	1	3.65
33	15.1821.01	Short type cylinders, PP, graduated, 500 ml	1	4.95
33	15.1822.01	Short type cylinders, PP, graduated, 1000 ml	1	27.50
33	15.1823.01	Short type cylinders, PP, graduated, 2000 ml	1	35.85
33	15.1332.02	Short type cylinders, TPX, graduated, 10 ml	1	6.35
33	15.1333.02	Short type cylinders, TPX, graduated, 25 ml	1	7.05
33	15.1334.01	Short type cylinders, TPX, graduated, 50 m	1	9.75
33	15.1335.01	Short type cylinders, TPX, graduated, 100 ml	1	12.70
33	15.1336.01	Short type cylinders, TPX, graduated, 250 ml	1	25.40
33	15.1337.01	Short type cylinders, TPX, graduated, 500 ml	1	30.80

Page	Item	Description	Pieces per unit	Price per Unit CHF
33	15.1339.01	Short type cylinders, TPX, graduated, 2000 ml	1	58.80
33	15.1801.02	BEAKER PP.cc. 25	1	1.35
33	15.1802.02	BEAKER PP.cc. 50	1	1.45
33	15.1803.02	BEAKER PP.cc. 100	1	1.70
33	15.1805.02	BEAKER PP.cc. 250	1	1.95
33	15.1806.02	BEAKER PP. cc. 500	1	3.30
33	15.1808.02	BEAKER PP.cc.1000	1	6.35
33	15.1810.02	BEAKER PP.cc.2000	1	8.80
33	15.1811.01	BEAKER PP.cc.3000	1	32.45
33	15.1812.01	BEAKER PP.cc.5000	1	50.20
33	15.1801.01	BEAKER PP.BLUE cc. 25	1	1.60
33	15.1802.01	BEAKER PP.BLUE cc. 50	1	1.70
33	15.1803.01	BEAKER PP.BLUE cc. 100	1	1.75
33	15.1805.01	BEAKER PP.BLUE cc. 250	1	2.45
33	15.1806.01	BEAKER PP.BLUE cc. 500	1	3.95
33	15.1808.01	BEAKER PP.BLUE cc.1000	1	5.65
33	15.1810.01	BEAKER PP.BLUE cc.2000	1	10.10
33	15.1542.01	BEAKER TPX.cc. 50	1	7.70
33	15.1543.01	BEAKER TPX. cc. 100	1	7.90
33	15.1545.01	BEAKER TPX.cc.250	1	9.15
33	15.1546.01	BEAKER TPX.cc. 500	1	18.80
33	15.1548.01	BEAKER TPX.cc 1000	1	23.75
33	15.1329.01	BEAKER TPX. cc.2000	1	33.15
33	15.1330.01	BEAKER TPX.cc. 3000	1	51.25
33	15.1331.02	BEAKER TPX. cc.. 5000	1	81.85
34	15.1460.01	Erlenmeyer of PP with screw cap, 50 ml	1	6.10
34	15.1461.01	Erlenmeyer of PP with screw cap, 100 ml	1	6.70
34	15.1463.01	Erlenmeyer of PP with screw cap, 250 ml	1	8.80
34	15.1464.02	Erlenmeyer of PP with screw cap, 500 ml	1	11.65
34	15.1465.01	Erlenmeyer of PP with screw cap, 1000 ml	1	14.60
34	15.1480.08	PP volumetric flask with HDPE stopper 25 ml	1	8.80
34	15.1481.07	PP volumetric flask with HDPE stopper 50 ml	1	11.30
34	15.1482.05	PP volumetric flask with HDPE stopper 100 ml	1	12.10
34	15.1483.01	PP volumetric flask with HDPE stopper 250 ml	1	14.70
34	15.1485.05	PP volumetric flask with HDPE stopper 500 ml	1	15.50
34	15.1486.01	PP volumetric flask with HDPE stopper 1000 ml	1	18.25
34	15.4801.09	Plastic stopper 10/19, HDPE	1	0.40
34	15.4811.01	Plastic stopper 12/21, HDPE	1	0.50
34	15.4831.01	Plastic stopper 14/23, HDPE	1	0.65
34	15.4851.01	Plastic stopper 19/26, HDPE	1	0.75
34	15.4861.01	Plastic stopper 24/29, HDPE	1	0.90
36	15.8110.01	PVC tubing mm 3X5	20 m	14.00
36	15.8110.02	PVC tubing mm 4X6	20 m	17.00
36	15.8110.03	PVC tubing mm 5X8	20 m	31.00
36	15.8110.04	PVC tubing mm 6X9	20 m	36.00
36	15.8110.05	PVC tubing mm 8X12	20 m	60.00
36	15.8110.06	PVC tubing mm 10X14	20 m	71.00
36	15.8110.07	PVC tubing mm 12X17	20 m	87.00
36	15.8120.01	Silicone tubing mm: 1X3	10 m	14.00
36	15.8120.02	Silicone tubing mm: 2X4	10 m	17.50
36	15.8120.03	Silicone tubing mm: 3X5	10 m	20.00
36	15.8120.04	Silicone tubing mm: 4X6	10 m	23.50
36	15.8120.05	Silicone tubing mm: 5X8	10 m	38.50
36	15.8120.06	Silicone tubing mm: 6X9	10 m	44.00
36	15.8120.07	Silicone tubing mm: 7X11	10 m	70.50
36	15.8120.08	Silicone tubing mm: 8X12	10 m	75.50
36	15.8120.09	Silicone tubing mm: 9X13	10 m	88.00
36	15.8120.10	Silicone tubing mm: 10X14	10 m	96.00
36	15.8121.03	Silicone tubing mm: 12X17	10 m	136.00
36	15.8121.04	Silicone tubing mm: 14X19	10 m	156.00
36	15.8130.01	PTFE tubing mm: 2X4	5 m	42.00
36	15.8130.02	PTFE tubing mm: 3X5	5 m	57.00
36	15.8130.03	PTFE tubing mm: 4X6	5 m	67.00
36	15.8130.04	PTFE tubing mm: 6X8	5 m	96.75
36	15.8130.05	PTFE tubing mm: 7X9	5 m	108.75
36	15.8130.06	PTFE tubing mm: 8X10	5 m	122.25
36	15.8130.07	PTFE tubing mm: 9X11	5 m	136.50
36	15.8130.08	PTFE tubing mm: 10X12	5 m	154.50
36	15.8130.09	PTFE tubing mm: 12X14	5 m	188.25
36	15.8140.01	Red rubber tubing mm: 3X5	10 m	37.50
36	15.8140.02	Red rubber tubing mm: 4X7	10 m	47.00
36	15.8140.03	Red rubber tubing mm: 5X8	10 m	44.50
36	15.8140.04	Red rubber tubing mm: 6X9	10 m	51.00
36	15.8140.05	Red rubber tubing mm: 6X12	10 m	104.50
36	15.8140.06	Red rubber tubing mm: 8X12	10 m	90.00
36	15.8140.07	Red rubber tubing mm: 9X13	10 m	89.50
36	15.8140.08	Red rubber tubing mm: 10X14	10 m	97.50
36	15.8140.09	Red rubber tubing mm: 12X17	10 m	148.50
36	15.8150.01	Latex tubing mm: 4X7	10 m	70.00

valid until 31.12.2017

Page	Item	Description	Pieces per unit	Price per Unit CHF
36	15.8150.03	Latex tubing mm: 8X12	10 m	121.00
36	15.8150.04	Latex tubing mm: 12X17	10 m	279.50
36	15.8160.01	Para tubing mm: 5X15	6 m	306.30
36	15.8160.02	Para tubing mm: 6X18	6	338.70
36	15.8160.03	Para tubing mm: 8X16	6	274.20
36	15.8160.04	Para tubing mm: 12X21	6	322.50
37	15.8071.01	Mono bulb with valve and tubes	1	25.05
37	15.8072.01	Vacuum bulb with two valves	1	60.20
37	15.8073.09	Richardson bulb for automatic burettes, Ø 46 mm	1	42.85
37	15.8073.10	Richardson bulb for automatic burettes, Ø 53 mm	1	44.55
37	15.8060.09	Latex gas sampling balloon mm 130	1	15.00
37	15.8060.10	Latex gas sampling balloon mm 150	1	19.50
37	15.8060.11	Latex gas sampling balloon mm 160	1	23.15
37	15.8060.13	Latex gas sampling balloon mm 200	1	26.80
37	15.8060.14	Latex gas sampling balloon mm 230	1	34.95
37	15.8060.15	Latex gas sampling balloon mm 250	1	43.15
37	15.8050.01	Reversible silicone stopper mm 7	30	51.00
37	15.8050.02	Reversible silicone stopper mm 10	30	57.00
37	15.8050.03	Reversible silicone stopper mm 13	30	97.50
37	15.8050.04	Reversible silicone stopper mm 15	20	84.00
37	15.8050.05	Reversible silicone stopper mm 16	20	92.00
37	15.8050.06	Reversible silicone stopper mm 20	10	64.50
37	15.8050.07	Reversible silicone stopper mm 24	10	108.00
39	15.8020.01	Silicone conical stopper mm 9	10	11.00
39	15.8020.02	Silicone conical stopper mm 12	10	13.00
39	15.8020.03	Silicone conical stopper mm 13	10	13.50
39	15.8020.04	Silicone conical stopper mm 14	10	14.00
39	15.8020.05	Silicone conical stopper mm 16	10	15.00
39	15.8020.06	Silicone conical stopper mm 18	10	16.00
39	15.8020.07	Silicone conical stopper mm 20	10	22.50
39	15.8020.08	Silicone conical stopper mm 23	10	33.50
39	15.8020.09	Silicone conical stopper mm 26	10	46.50
39	15.8020.10	Silicone conical stopper mm 30	10	65.00
39	15.8021.06	Silicone conical stopper mm 32	5	54.00
39	15.8021.07	Silicone conical stopper mm 37	5	82.25
39	15.8021.08	Silicone conical stopper mm 42	5	107.50
39	15.8021.09	Silicone conical stopper mm 45	5	119.50
39	15.8021.10	Silicone conical stopper mm 48	5	154.50
39	15.8021.11	Silicone conical stopper mm 54	5	170.75
39	15.8021.12	Silicone conical stopper mm 60	5	219.25
39	15.8030.01	Conical red rubber stopper, mm 9	25	5.00
39	15.8030.02	Conical red rubber stopper, mm 12	25	6.25
39	15.8030.03	Conical red rubber stopper, mm 13	25	7.50
39	15.8030.04	Conical red rubber stopper, mm 15	25	7.50
39	15.8030.05	Conical red rubber stopper, mm 16	25	8.75
39	15.8030.06	Conical red rubber stopper, mm 18	25	10.00
39	15.8030.07	Conical red rubber stopper, mm 20	25	11.25
39	15.8030.08	Conical red rubber stopper, mm 23	25	21.25
39	15.8030.09	Conical red rubber stopper, mm 26	25	28.75
39	15.8030.10	Conical red rubber stopper, mm 30	25	33.75
39	15.8031.01	Conical red rubber stopper, mm 32	25	51.25
39	15.8031.02	Conical red rubber stopper, mm 37	5	19.50
39	15.8031.03	Conical red rubber stopper, mm 42	5	23.50
39	15.8031.04	Conical red rubber stopper, mm 45	5	27.00
39	15.8031.05	Conical red rubber stopper, mm 48	5	33.00
39	15.8031.06	Conical red rubber stopper, mm 54	5	40.00
39	15.8031.07	Conical red rubber stopper, mm 60	5	57.75
39	15.8031.08	Conical red rubber stopper, mm 65	5	89.50
39	15.8040.03	Conical red rubber stopper with 1 hole, mm 12	25	7.50
39	15.8040.04	Conical red rubber stopper with 1 hole, mm 13	25	7.50
39	15.8040.05	Conical red rubber stopper with 1 hole, mm 14	25	8.75
39	15.8040.06	Conical red rubber stopper with 1 hole, mm 16	25	8.75
39	15.8040.07	Conical red rubber stopper with 1 hole, mm 18	25	12.50
39	15.8040.08	Conical red rubber stopper with 1 hole, mm 20	25	13.75
39	15.8040.09	Conical red rubber stopper with 1 hole, mm 23	25	22.50
39	15.8040.10	Conical red rubber stopper with 1 hole, mm 26	25	33.75
39	15.8040.11	Conical red rubber stopper with 1 hole, mm 30	25	38.75
39	15.8040.12	Conical red rubber stopper with 1 hole, mm 32	25	53.75
39	15.8041.05	Conical red rubber stopper with 1 hole, mm 37	5	19.75
39	15.8041.06	Conical red rubber stopper with 1 hole, mm 42	5	25.00
39	15.8041.07	Conical red rubber stopper with 1 hole, mm 45	5	28.00
39	15.8041.08	Conical red rubber stopper with 1 hole, mm 48	5	32.25
39	15.8041.09	Conical red rubber stopper with 1 hole, mm 54	5	38.75
39	15.8041.10	Conical red rubber stopper with 1 hole, mm 60	5	71.25
39	15.8041.11	Conical red rubber stopper with 1 hole, mm 65	5	92.50
39	15.8045.06	Conical red rubber stopper with 2 holes, mm 18	25	20.00
39	15.8045.07	Conical red rubber stopper with 2 holes, mm 20	25	15.00
39	15.8045.08	Conical red rubber stopper with 2 holes, mm 23	25	22.50
39	15.8045.09	Conical red rubber stopper with 2 holes, mm 26	25	31.25

Page	Item	Description	Pieces per unit	Price per Unit CHF
39	15.8045.11	Conical red rubber stopper with 2 holes, mm 32	25	51.25
39	15.8045.12	Conical red rubber stopper with 2 holes, mm 37	5	20.75
39	15.8045.13	Conical red rubber stopper with 2 holes, mm 42	5	25.75
39	15.8045.14	Conical red rubber stopper with 2 holes, mm 45	5	30.25
39	15.8045.15	Conical red rubber stopper with 2 holes, mm 48	5	34.50
39	15.8046.10	Conical red rubber stopper with 2 holes, mm 54	5	40.00
39	15.8046.11	Conical red rubber stopper with 2 holes, mm 60	5	75.50
39	15.8046.12	Conical red rubber stopper with 2 holes, mm 65	5	91.25
39	15.8010.01	PTFE magnetic stirring bar mm 6 X 10	10	19.00
39	15.8010.02	PTFE magnetic stirring bar mm 6 X 15	10	23.00
39	15.8010.03	PTFE magnetic stirring bar mm 6 X 20	10	26.50
39	15.8010.04	PTFE magnetic stirring bar mm 6 X 25	10	32.50
39	15.8010.05	PTFE magnetic stirring bar mm 6 X 30	10	36.00
39	15.8010.06	PTFE magnetic stirring bar mm 6 X 35	10	45.00
39	15.8010.07	PTFE magnetic stirring bar mm 8 X 40	10	54.00
39	15.8010.08	PTFE magnetic stirring bar mm 8 X 45	10	68.00
39	15.8010.09	PTFE magnetic stirring bar mm 8 X 50	10	73.00
39	15.8010.10	PTFE magnetic stirring bar mm 9 x 60	10	82.50
39	15.8011.01	PTFE magnetic stirring bar mm 9 x 70	10	115.00
39	15.8075.01	Latex bulbs for glass pasteur pipettes	100	70.00
40	15.3001.01	10 Lt stack crate, white, HDPE	1	16.15
40	15.3003.01	16 Lt stack crate, white, HDPE	1	27.70
40	15.3005.02	16 Lt stack crate, white, HDPE, with holes	1	25.55
40	15.3004.02	55 Lt stack crate, white, HDPE	1	50.05
40	15.1602.06	20 Lt stack crate, white, HDPE	1	21.15
41	15.1603.01	50 Lt stack crate, white, HDPE	1	40.65
41	15.1604.06	80 Lt stack crate, white, HDPE	1	58.15
41	15.1601.09	100 Lt stack crate, white, HDPE	1	111.90
42	15.1780.01	20 Lt cool box, blue, HDPE	1	68.00
42	15.1782.01	32 Lt cool box, blue, HDPE	1	91.95
42	15.1783.01	45 Lt cool box, blue, HDPE	1	154.00
42	15.1066.01	PP scoop, 2000 ml	60	13.10
42	15.2147.01	Large cup, white, HDPE, 3000 ml	1	13.45
42	15.1062.02	Round form scoop, 18 cm	1	2.60
42	15.1063.02	Round form scoop, 24 cm	1	2.85
42	15.1064.03	Round form scoop, 30 cm	1	3.20
42	15.1065.01	Round form scoop, 36 cm	1	4.40
42	15.1067.02	Round form scoop, 40 cm	1	5.85
43	15.2032.01	Spatulas in HDPE, white, 600 mm	1	82.30
43	15.2033.01	Spatulas in HDPE, white, 800 mm	1	89.05
43	15.2035.01	Spatulas in HDPE, white, 1000 mm	1	104.15
43	15.2038.01	Spatulas in HDPE, white, 1200 mm	1	110.90
43	15.2039.01	Spatulas in HDPE, white, 1400 mm	1	124.30
43	15.1955.01	Suction pump for acids en PP, 250 ml	1	329.30
43	15.1935.01	Pump HDPE/PVC	1	20.50
44	15.1400.08	Male bed bottle, HDPE	1	7.95
44	15.1401.06	PP bed basin	1	13.45
44	15.1402.03	Bottle carrier, red, HDPE, for 6 bottles of max. Ø 95 mm	1	8.85



committed to

HUBERLAB.

committed to science

HUBERLAB. AG
Industriestrasse 123
4147 Aesch

T 061 717 99 77
F 061 711 93 42

info@huberlab.ch
www.huberlab.ch

Prices net, excl. VAT / free delivery from CHF 300.00 net