

Nephelo Culture Flask Wheaton

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Consumables](#) / [Containers](#) / [Flasks](#) / Nephelo Culture Flask Wheaton

Product Data

- For use with turbidity meter or nephelometer
- Scratch and striation free side tube
- Available with or without cleanout port
- Autoclavable borosilicate glass that conforms to USP Type I and ASTM E 438 Type I, Class A requirements
- 38-430 top cap and 20-415 sidearm cap have white styrene butadiene rubber liner

Item No.	volume	vol.	description	execution	Price
8.3510.54	500 ml	350 ml	Sidearm Size: 14 x 130 mm	without Cleanout Port	CHF 86.20
8.3514.54	300 ml	200 ml	Sidearm Size: 14 x 140 mm	with Cleanout Port	CHF 204.90
8.3514.56	300 ml	200 ml	Sidearm Size: 19 x 140 mm	with Cleanout Port	CHF 175.85
8.3514.86	500 ml	250 ml	Sidearm Size: 19 x 140 mm	with Cleanout Port	CHF 218.10

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 02.06.2024 at 01:13 oclock.

www.huberlab.ch