

# HYDRANAL™ - Water Standard Oil, Oil-based standard for Karl Fischer titration (water

**Honeywell  
Fluka™**

Produkt im eShop

Startseite / Covid-19 / Covid-19 Honeywell / Covid-19 Honeywell - Quality control for drugs and vaccines / Covid-19 Honeywell - Quality control for drugs and vaccines - Standards / HYDRANAL™ - Water Standard Oil, Oil-based standard for Karl Fischer titration (water

## Produktdaten

CAS-Nr	8042-47-5
EG-Nr	232-455-8
Lagerung	Room Temp
LGK	10
Sicherheitshinweis(e)	P280
WGK	1

Art.-Nr.	Inhalt	Verpackung	Preis
34694-80ML	80 Milliliter	80ML Ampoule	CHF 219.06

\* Die Preise sind freibleibend und verstehen sich als Abgabepreise in Schweizer Franken ohne Mehrwertsteuer (MWST), sowie exklusive aller sonstigen Gebühren, Abgaben und Steuern. Die Preise, die im eShop angezeigt werden, sind aufgrund der regelmässigen Aktualisierung möglicherweise aktueller als auf der PDF-Datei.

\*\* Bitte beachten Sie, dass bei der Bestellung von Chemikalien und Reinigungsmittel Transport - und Verpackungskosten für Gefahrgüter sowie gesetzlich vorgeschriebene Gebühren berechnet werden. Diese werden Ihnen auf der Auftragsbestätigung, die Sie ergänzend zur Bestellbestätigung erhalten, detailliert ausgewiesen.

\*\*\* Weitere Informationen wie Technische Informationen und Sicherheitsdatenblätter finden Sie online in unserem eShop

\*\*\*\* Die PDF-Datei wurde erzeugt auf [www.huberlab.ch](http://www.huberlab.ch) am 15.05.2024 um 22:02 Uhr.