

# AQUAMETRIC Solvent for volumetric analysis



Product in eShop

Home / Assortment / Chemicals / Laboratory Chemicals / Mass solutions and standards / AQUAMETRIC Solvent for volumetric analysis

## Product Data

HS-No	38220000
LGK	3
Hazard Statement(s)	H225-H301-H311-H314-H331-H370
Safety Instructions	P210-P233-P240-P241-P242-P243-P260-P261-P264-P270-P271-P280-P301+P310-P301+P330+P331-P302+P352-P303+P361+P353-P304+P340-P305+P351+P338-P307+P311-P310-P311-P312-P321-P322-P330-P338-P361-P363-P370+P378-P403+P233-P403+P235-P405-P501
Signal Word	Gefahr
Class / PG	3(6.1,8)/II
UN-No.	UN3286
WHC	1



Item No.	Content	Packaging	Price
285817.1611	1'000 milliliter	Glasflasche	CHF 112.90
285817.1612	2.5 liter	Glas	CHF 235.40

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 03.05.2024 at 07:31 oclock.