

Sigma 2-16P bench-top centrifuge



Product in eShop

[Home](#) / [Products](#) / Sigma 2-16P bench-top centrifuge

Product Data

Item No.	Description	UNSPSC	ECLASS	Price
17.0101.07	Sigma 2-7 Laboratory centrifuge	41100000	32010000	on demand
17.0101.55	Sigma 2-16P Labortischzentrifuge	41100000	32010000	on demand
17.0101.71	Swingout-rotor 4 x 100 ml w/o buckets	41100000	32010000	on demand
17.0101.72	Round bucket 100ml with adapter for roundbottom-tubes	41100000	32010000	on demand
17.0101.73	Rubber cushion for 50&100ml	41100000	32010000	on demand
17.0101.78	Glass tube with round bottom, 100 ml w/o buckets	41100000	32010000	on demand

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 14:11 o'clock.