

Ice bucket PUR with lid, red, LABSOLUTE®;

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Consumables](#) / [Containers](#) / [Ice tubs & cool boxes](#)
/ Ice bucket PUR with lid, red, LABSOLUTE®;

Product Data

The red LABSOLUTE® ice buckets are made of polyurethane foam and ideal for the cooling of samples, solutions or cultures. The ice buckets have excellent insulating properties and are suitable for use with dry ice, water ice, brine and acetone.

- Extremely robust and light
- With lid

Item No.	vol. l	dimensions (l x w x h)	colour	Price
32.2137.70	2.5 l	283 x 283 x 130 mm	red	CHF 112.27
32.2137.80	4.5 l	280 x 280 x 185 mm	red	CHF 192.02

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 28.05.2024 at 20:06 oclock.