


Test tube rack, wire/nylon LLG-Labware

[Product in eShop](#)

Home / [Labortops Q2 2026](#) / Test tube rack, wire/nylon LLG-Labware

Product Data

Made of wire with white nylon coating. Suitable for safe storage and transport of test tubes and centrifuge tubes. The second tier provides an optimal hold. Opening size 18 x 18 mm (4.711 546 = 30 x 30 mm).

- Colour: White
- Material: Wire

Item No.	Amount/Quantity	Ø opening	dimensions (w x d x h)	Price
600.7115.42	12 positions	18 mm	54 x 130 x 70 mm	CHF 8.50 CHF 13.40
600.7115.44	48 positions	18 mm	94 x 248 x 70 mm	CHF 12.05 CHF 19.00
600.7115.49	100 positions	18 mm	208 x 215 x 70 mm	CHF 28.35 CHF 44.65

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 21:37 oclock.

www.huberlab.ch