

# ICP Multi Element Tuning Standards Reagecon

[Product in eShop](#)

[Home](#) / [Anion&Cation](#) / ICP Multi Element Tuning Standards Reagecon

## Product Data

Reagecon have been manufacturing Inorganic Standards, Controls and Calibrators for Spectroscopy for almost two decades. During that time, the company has established itself as the most reliable primary producer of top quality standards. The products are manufactured, tested and stabilised to such a high level, that they can be used on all of the instruments. Reagecon's Multi Element Standards for ICP, ICP-MS are manufactured from ultra-pure metals or salts. Raw materials used for these standards are assayed by titration and ICP-MS prior to manufacture. Separate CRM's are used to control or calibrate the titration and ICP-MS respectively. Products are manufactured gravimetrically using the mass balance approach: 100% less the sum of all impurities (w/w).

Item No.	Content	Packaging	Price
500.REICPCAL4X5	500 milliliter	HDPE Triple-leached	CHF 599.02
500.REICPTUNE10B	100 milliliter	HDPE Triple-leached	CHF 265.15
500.REICPTUNE11A	100 milliliter	HDPE Triple-leached	CHF 562.55
500.REICPTUNE2	100 milliliter	HDPE Triple-leached	CHF 380.02
500.REICPTUNE22B	100 milliliter	HDPE Triple-leached	CHF 946.00
500.REICPTUNE23A	100 milliliter	HDPE Triple-leached	CHF 942.22
500.REICPTUNE25A	100 milliliter	HDPE Triple-leached	CHF 1'112.50
500.REICPTUNE2A	100 milliliter	HDPE Triple-leached	CHF 276.25
500.REICPTUNE33A	100 milliliter	HDPE Triple-leached	CHF 930.85
500.REICPTUNE3A	100 milliliter	HDPE Triple-leached	CHF 336.79
500.REICPTUNE4C	100 milliliter	HDPE Triple-leached	CHF 240.75
500.REICPTUNE5	100 milliliter	HDPE Triple-leached	CHF 700.08
500.REICPTUNE5A	100 milliliter	HDPE Triple-leached	CHF 374.60
500.REICPTUNE5C	100 milliliter	HDPE Triple-leached	CHF 394.50
500.REICPTUNE6	100 milliliter	HDPE Triple-leached	CHF 438.94
500.REICPTUNE6A	100 milliliter	HDPE Triple-leached	CHF 404.91
500.REICPTUNE7C	100 milliliter	HDPE Triple-leached	CHF 350.75
500.REICPTUNE7E	100 milliliter	HDPE Triple-leached	CHF 205.40
500.REICPTUNE7F	100 milliliter	HDPE Triple-leached	CHF 146.54
500.REICPTUNE9A	100 milliliter	HDPE Triple-leached	CHF 503.25
500.REICPTUNE9B	100 milliliter	HDPE Triple-leached	CHF 503.25

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 21.04.2026 at 14:11 oclock.

[www.huberlab.ch](http://www.huberlab.ch)