

# Colour Standards Reagecon

[Product in eShop](#)[Home](#) / [Brands](#) / [Reagecon Diagnostics](#) / [Colour Standards Reagecon](#)

## Product Data

Reagecon offers the full range of ASTM, Saybolt, Platinum-Cobalt (Hazen), Gardner, European Pharmacopeia and United States Pharmacopeia Colour Standards for use with ASTM, APHA, ACS, Europeann and United States Pharmacopeia standard methods.

Item No.	Content	Packaging	Price
500.ASTMA051	100 milliliter	Amber Glass Winchester Bottle	CHF 101.70
500.ASTMA055	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 597.05
500.ASTMA101	100 milliliter	Amber Glass Winchester Bottle	CHF 101.70
500.ASTMA105	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 322.75
500.ASTMA301	100 milliliter	Amber Glass Winchester Bottle	CHF 101.70
500.ASTMA305	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 322.75
500.ASTMA501	100 milliliter	Amber Glass Winchester Bottle	CHF 132.25
500.ASTMA505	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 597.05
500.ASTMA701	100 milliliter	Amber Glass Winchester Bottle	CHF 161.00
500.ASTMA705	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 597.05
500.GARD011	100 milliliter	Amber Glass Winchester Bottle	CHF 271.20
500.GARD02	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 237.10
500.GARD021	100 milliliter	Amber Glass Winchester Bottle	CHF 188.15
500.GARD031	100 milliliter	Amber Glass Winchester Bottle	CHF 146.60
500.GARD04	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 237.10
500.GARD041	100 milliliter	Amber Glass Winchester Bottle	CHF 188.15
500.GARD051	100 milliliter	Amber Glass Winchester Bottle	CHF 146.60
500.GARD06	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 237.10

500.GARD061	100 milliliter	Amber Glass Winchester Bottle	CHF 188.15
500.GARD071	100 milliliter	Amber Glass Winchester Bottle	CHF 146.60
500.GARD08	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 237.10
500.GARD081	100 milliliter	Amber Glass Winchester Bottle	CHF 188.15
500.GARD091	100 milliliter	Amber Glass Winchester Bottle	CHF 146.60
500.GARD10	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 237.10
500.GARD101	100 milliliter	Amber Glass Winchester Bottle	CHF 101.70
500.GARD111	100 milliliter	Amber Glass Winchester Bottle	CHF 146.60
500.GARD12	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 237.10
500.GARD121	100 milliliter	Amber Glass Winchester Bottle	CHF 188.15
500.GARD131	100 milliliter	Amber Glass Winchester Bottle	CHF 271.20
500.GARD14	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 237.10
500.GARD141	100 milliliter	Amber Glass Winchester Bottle	CHF 188.15
500.GARD151	100 milliliter	Amber Glass Winchester Bottle	CHF 146.60
500.GARD16	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 294.85
500.GARD161	100 milliliter	Amber Glass Winchester Bottle	CHF 188.15
500.GARD171	100 milliliter	Amber Glass Winchester Bottle	CHF 146.60
500.GARD181	100 milliliter	Amber Glass Winchester Bottle	CHF 146.60
500.HAZ0	1 liter	HDPE Bottle GL45	CHF 124.10
500.HAZ001	100 milliliter	HDPE Bottle	CHF 63.65
500.HAZ10	1 liter	HDPE Bottle GL45	CHF 124.10
500.HAZ100	1 liter	HDPE Bottle GL45	CHF 275.85
500.HAZ1000	1 liter	HDPE Bottle GL45	CHF 1'091.90
500.HAZ1001	100 milliliter	HDPE Bottle	CHF 63.20
500.HAZ15	1 liter	HDPE Bottle GL45	CHF 229.65
500.HAZ150	1 liter	HDPE Bottle GL45	CHF 263.85
500.HAZ20	1 liter	HDPE Bottle GL45	CHF 229.65
500.HAZ200	1 liter	HDPE Bottle GL45	CHF 502.05
500.HAZ2005	5 liter	HDPE Container	CHF 975.87
500.HAZ25	1 liter	HDPE Bottle GL45	CHF 124.10
500.HAZ250	1 liter	HDPE Bottle GL45	CHF 275.85
500.HAZ30	1 liter	HDPE Bottle GL45	CHF 408.25
500.HAZ40	1 liter	HDPE Bottle GL45	CHF 220.70
500.HAZ5	1 liter	HDPE Bottle GL45	CHF 124.15
500.HAZ50	1 liter	HDPE Bottle GL45	CHF 220.70
500.HAZ500	1 liter	HDPE Bottle GL45	CHF 468.95
500.HAZ500-500ml	500 milliliter	HDPE Bottle GL45	CHF 250.30
500.HAZ50001	100 milliliter	HDPE Bottle	CHF 131.30
500.HAZ5001	100 milliliter	HDPE Bottle	CHF 101.65
500.HAZ5005	5 liter	HDPE Jerrycan	CHF 1'031.50
500.HAZ501	100 milliliter	HDPE Bottle	CHF 61.00
500.HAZ505	5 liter	HDPE Container	CHF 508.62
500.HAZ80	1 liter	HDPE Bottle GL45	CHF 510.35
500.HAZ8001	100 milliliter	HDPE Bottle	CHF 75.55

500.SAYN105	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 232.60
500.SAYN151	100 milliliter	Amber Glass Winchester Bottle	CHF 184.60
500.SAYN155	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 232.60
500.SAYP01	100 milliliter	Amber Glass Winchester Bottle	CHF 184.60
500.SAYP05	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 240.00
500.SAYP121	100 milliliter	Amber Glass Winchester Bottle	CHF 184.60
500.SAYP125	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 240.00
500.SAYP155	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 240.00
500.SAYP191	100 milliliter	Amber Glass Winchester Bottle	CHF 184.60
500.SAYP195	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 240.00
500.SAYP251	100 milliliter	Amber Glass Winchester Bottle	CHF 184.60
500.SAYP255	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 232.60
500.SAYP301	100 milliliter	Amber Glass Winchester Bottle	CHF 99.75
500.SAYP305	500 milliliter	Amber Glass Winchester Bottle GL45	CHF 232.60

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 22.04.2026 at 07:03 o'clock.