

## Urine bottle for women

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Medicine & Healthcare](#) / [Sampling](#) / [Urine](#) / [Urine collection](#) / [Urine bottle for women](#)

### Product Data

1 litre capacity, translucent, graduated, total length approx. 30 cm, heat-resistant up to 90°C, autoclavable up to 130°C, may be disinfected, suitable for all common urine bottle holders, with attached cap made of white polyethylene

Item No.	Price
3.4104.02	CHF 9.25

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 01.05.2024 at 19:36 o'clock.