

# Laboratory Bottle DURAN® YOUTILITY DWK

[Product in eShop](#)

Home / [Labortops Q2 2026](#) / Laboratory Bottle DURAN® YOUTILITY DWK

## Product Data

### DURAN® YOUTILITY Bottle System

The DURAN Group introduces the new bottle system for the 21st century laboratory. The new YOUTILITY bottle system was developed by the DURAN Group, who originated the screw threaded GL 45 DURAN® borosilicate glass laboratory bottle in 1972. Since then the DURAN® GL 45 Laboratory Bottle has become the bottle of choice for laboratories around the world.

Now, the DURAN® YOUTILITY bottle represents a new generation of glass laboratory bottles for use in a wide and diverse range of applications. The innovation does not stop at the improved ergonomic bottle shape, but extends to dedicated accessories that incorporate many new innovative features to improve handling, sample identification, and ease of use. The new YOUTILITY bottles and accessories help to make laboratory work easier, safer, more economical and fun.

### DURAN® YOUTILITY Laboratory Bottle GL 45

- Specially shaped gripping zones on both sides of the bottle for easier and safer handling.
- With the new YOUTILITY bottle thread opening or closing the bottle is significantly faster. The thread is fully compatible with DIN GL 45 closures and other accessories.
- The slimmer YOUTILITY bottle shape allows a more optimal use of limited space in autoclaves and laboratory refrigerators.
- A pre-defined labelling area is compatible with the dedicated DURAN® self-adhesive YOUTILITY bottle labels.
- Nominal volume is shown at the top of the easy-to-read graduation scale for fast determination of the volumes.
- Glass bottle body is moulded from the tried and tested DURAN® borosilicate 3.3 pharmacopoeial Type 1 neutral glass. DURAN® glass offers very good chemical resistance and high temperature resistance.
- The YOUTILITY bottle has an integrated Retrace Code to allow full product traceability. The lot specific certificate of conformance is available at: [www.duran-group.com](http://www.duran-group.com).

- Each DURAN® YOUTILITY bottle is supplied as a complete system, with a pouring ring (cyan polypropylene) and a GL 45 cap (cyan polypropylene).

Item No.	volume	d x h	description	Price
8.5817.28	125 ml	55 x 124 mm	Youtility bottle	CHF 24.85 <del>CHF 49.80</del>
8.5817.36	250 ml	66 x 158 mm	Youtility bottle	CHF 27.70 <del>CHF 55.60</del>
8.5817.44	500 ml	78 x 193 mm	Youtility bottle	CHF 34.50 <del>CHF 69.20</del>
8.5817.54	1000 ml	93 x 253 mm	Youtility bottle	CHF 48.70 <del>CHF 97.80</del>
8.5827.28	125 ml	55 x 124 mm	Youtility Laborflasche Braunglas	CHF 81.00 <del>CHF 162.60</del>
8.5827.36	250 ml	66 x 158 mm	Youtility Laborflasche Braunglas	CHF 112.10 <del>CHF 224.80</del>
8.5827.44	500 ml	78 x 193 mm	Youtility Laborflasche Braunglas	CHF 144.90 <del>CHF 290.60</del>
8.5827.54	1000 ml	93 x 253 mm	Youtility Laborflasche Braunglas	CHF 177.00 <del>CHF 356.40</del>

## Accessories

Item No.	Description	Price
8.5817.03		from CHF 41.60
8.5817.01	Accessories for Laboratory Bottle DURAN YOUTILITY DWK	from CHF 14.40
8.5817.14	DURAN®, YOUTILITY Screw Cap GL45	from CHF 36.00

## More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 21.04.2026 at 06:36 oclock.

[www.huberlab.ch](http://www.huberlab.ch)