



## Microscope slide boxes



Product in eShop

Home / [Assortment](#) / [Medicine & Healthcare](#) / [Microscopy](#) / [Slides](#) / Microscope slide boxes

### Product Data

Constructed of durable polypropylene, slide boxes hold either 25 or 100 standard slides (76x 26 mm)

- Hinged lids are held securely closed by latch lock
- Stacking feet maximize stability
- Slide Identification card
- No liner

Item No.	for number of slides	colour	Dimensions	Price
15.1599.55	25	blue	141 x 92 x 3 mm	CHF 7.50
15.1599.56	25	red	141 x 92 x 3 mm	CHF 7.50
15.1599.57	25	white	141 x 92 x 3 mm	CHF 7.50
15.1599.58	100	blue	169 x 210 x 37 mm	CHF 10.15
15.1599.59	100	red	169 x 210 x 37 mm	CHF 10.15
15.1599.60	100	white	169 x 210 x 37 mm	CHF 10.15

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 18.04.2026 at 20:31 o'clock.