


Biopsy cassettes in QuickLoad™ taped stacks, Unisette™

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Microscopie](#) / Biopsy cassettes in QuickLoad™ taped stacks, Unisette™

Product Data

Made of acetal polymer, molded from a special high density polymer, these cassettes keep specimens safely submerged in liquid and are totally resistant to the chemical action of histological solvents. The efficient flow-through slots maximise fluid exchange and ensure proper drainage. The one piece integral lid eliminates the need for separate steel lids. Lids can be opened and closed as often as necessary and they always relock securely without danger of specimen loss. Anterior writing area is at a 35° angle. M506 series is similar to M505 but specially designed to hold biopsy specimens during the embedding process.

- The snaplatch and hinge-lock design prevents early separation of base and lid and allows for one-handed operation
- Large labelling areas are located on two sides of the cassettes

Item No.	colour	Price
10.9506.02	white	CHF 183.00
10.9506.03	pink	CHF 106.48
10.9506.04	green	CHF 106.48
10.9506.05	yellow	CHF 106.48
10.9506.06	blue	CHF 212.71
10.9506.07	peach	CHF 106.48
10.9506.08	beige	CHF 106.48
10.9506.09	gray	CHF 106.48
10.9506.10	purple	CHF 106.48
10.9506.11	orange	CHF 106.48
10.9506.12	aqua	CHF 106.48

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 06:36 o'clock.

www.huberlab.ch