

Paraffin Block Mailer



Product in eShop



Home / [Assortment](#) / [Research & Development](#) / [Microscopie](#) / Paraffin Block Mailer

Product Data

Transport container made of PVC can be used for handling and shipping up to 6 blocks as needed, while having to keep only one model of block mailer in inventory.

- Secure and easy to use
- Transparent for easy viewing of contents
- Suitable for all regular models of tissue and biopsy cassettes
- Dimensions: 134 x 141 x 29 mm

Item No.	Price
15.9477.99	CHF 97.27

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 14:11 o'clock.