

Burettes, lateral stopcock SILBERBRAND, class B



Product in eShop

Home / [Assortment](#) / [Research & Development](#) / [Titration](#) / Burettes, lateral stopcock SILBERBRAND, class B

Product Data

Boro 3.3. To DIN EN ISO 385. Calibrated to deliver (TD, Ex). Schellbach stripe. Interchangeable needle-valve and glass key. Length approx. 800 mm.

Item No.	volume	grad.	execution	Price
1.1631.66	25 ml	0,05 ml	stopcock with PTFE spindle	CHF 139.15
1.1631.73	50 ml	0,1 ml	stopcock with glass key	CHF 57.75
1.1631.78	10 ml	0,02 ml	stopcock with PTFE spindle	CHF 139.15
1.1631.79	50 ml	0,1 ml	stopcock with PTFE spindle	CHF 139.15

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 29.04.2024 at 14:21 oclock.