

Reagent bottles opaque wide mouth

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Quality Assurance](#) / [Sampling](#) / [Narrow neck bottles](#) / [Reagent bottles opaque wide mouth](#)

Product Data

Made from opaque, pigmented PP for light-sensitive substances. According to DIN 12039. With NS stopper made of PP.

- Good chemical resistance
- Ideal for long-term storage of liquids
- Wider mouth for simple filling, also ideal for powders and paste-like materials
- Autoclavable at 121 °C (2 bar) according to DIN EN 285

Item No.	vol.	NS	Price
15.1017.08	500 ml	45 / 40	CHF 31.30

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 06:36 oclock.

www.huberlab.ch