

Erlenmeyer flasks wide mouth with screw caps Nalgene®



Product in eShop

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Consumables](#) / [Containers](#) / [Flasks](#) / [Erlenmeyer flasks](#) / Erlenmeyer flasks wide mouth with screw caps Nalgene®

Product Data

PMP (TPX), transparent with PP screw cap, ideal for titrating and for preparing and storing nutrient media

- Good chemical resistance
- Heat resistant up to +150 °C (up to +170 °C for short periods), autoclavable
- Good microwave permeability

Item No.	vol.	thread	Price
15.4109.45	125 ml	33	CHF 126.96
15.4109.46	250 ml	38	CHF 88.05
15.4109.47	500 ml	43	CHF 527.56
15.4109.48	1000 ml	53	CHF 190.02

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 16.05.2024 at 05:46 oclock.