

Microplate Centrifuge Greiner Bio-One

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Centrifugation](#) / [Centrifuges](#) / [Mini centrifuges](#) / [Microplate Centrifuge Greiner Bio-One](#)

Product Data

Two plates can be inserted confidently without sealing tape or caps. Upon closing the lid of the centrifuge, the rotor accelerates and centrifugal force 'swings' the plates into a horizontal position. With a g-force power of 600 x g, most samples can be spun down in 20 seconds or less. Simply open the lid and the automatic electric brake brings the rotor to a smooth stop in about 7 seconds. Plates can then easily be removed.

- Quickly spins down droplets in microplates
- Unique rotor design holds 2 plates and prevents spillage
- 50 % smaller than traditional centrifuges
- Compatible with most Greiner Bio-One microplates
- Speed 2550 rpm/600 x g, capacity 2 x PCR or microtiter plates

Item No.	description	Price
7.846 070	Microplate centrifuge, speed 2100 rpm/400 x g, capacity 2 x PCR or microtiter plates	CHF 1'838.99

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 04.05.2024 at 19:48 oclock.