

Folding Thermometer TLC 700 ebro

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Food technology & cuisine](#) / [Cuisine](#) / [Temperature measurement & data recording](#) / [Thermometers](#) / Folding Thermometer TLC 700 ebro

Product Data

- Small size easily fits in a pocket
- Foldable penetration probe for safe storage
- Waterproof (IP 65)
- According to EN 13485 (ability: S, T, E)
- Applications: For temperature measurement during transport, storage, and incoming goods inspection.

Technical data

- Measuring range: -30 °C...+220 °C
- Accuracy: ± 0.5 °C to -30 °C...+100 °C ± 1.0 °C for the remaining measurement range
- Resolution: 0.1 °C
- Operating Temperature: 0°...+50 °C
- Measured variables: Temperature
- IP Rating: IP 65
- Battery: Lithium button cell 3 V, Typ CR 2032
- Dimensions: 118 x 33 x 15 mm
- Weight approx.: 41,2 g
- Certificate: Factory calibration certificate

Item No.	temp.	division	Price
4.4332.10	-30 ... +220 °C	0,1 °C	CHF 67.00

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 11.05.2024 at 17:38 oclock.

www.huberlab.ch