

Desiccator SICCO Bohlender



Product in eShop



[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Vacuum Technology](#) / [Exsiccators](#) / [Exsiccators](#) / [Desiccator SICCO Bohlender](#)

Product Data

SICCO Desiccators are designed for storing or drying humidity sensitive products using silica gel. The tight fitting door protects the contents from contamination from the atmosphere. The controlled environment inside the desiccators are ideal for storing reference materials, retained samples, metallographic specimen, tobacco and DNA samples.

- Aluminium frame with panels made of acrylic glass
(for Star-Vitrum made of Borosilicate glass 3.3)
- Door with magnetic catch and circular rubber seal
- Non-slip rubber feet
- Easy to read electronic digital hygrometer with temperature and humidity display
- Variable height shelves made of acrylic glass
(for Star-Vitrum, Star-Big and Star-Auto-Big made of stainless steel)
- Telescopic shelves
- Desiccant tray
- Silica gel included
- Raster and rails made metal-free PA/fiberglass
- Rails with level numbering

* per chamber

Item No.	type no.	inner dimensions	vol. l	heat resistance	Price
15.1877.07	Auto-Star-Exsikkator	260 x 480 x 330 mm	51 l	-20 - +70 °C °C	CHF 2'470.00
15.1879.07	Star-Desiccator Protect	260 x 480 x 330 mm	51 l	-20 - +70 °C	CHF 1'719.00
15.1880.01	Mini Premium (1 compartment)	208 x 159 x 180 mm	6 l	-20 - +70 °C	CHF 249.10
15.1880.02	Mini Premium (2 compartments)	* 208 x 159 x 180 mm	* 6 l	-20 - +70 °C	CHF 423.00

15.1880.03	Mini Premium (3 compartments)	* 208 x 159 x 180 mm	* 6 l	-20 - +70 °C	CHF 629.00
15.1880.10	Star	260 x 480 x 330 mm	51 l	-20 - +70 °C	CHF 1'380.00
15.1880.11	Star-Vitrum	260 x 480 x 330 mm	51 l	-20 - +150 °C	CHF 1'945.00
15.1880.12	Big-Star	495 x 500 x 540 mm	156 l	-20 - +70 °C	CHF 2'877.00
15.1880.13	Big-Star Vitrum	495 x 500 x 540 mm	156 l	-20 - +150 °C	CHF 4'032.00
15.1880.14	Maxi 1	495 x 1030 x 540 mm	311 l	-20 - +70 °C	CHF 4'766.00
15.1880.15	Maxi 1 Vitrum	495 x 1030 x 540 mm	311 l	-20 - +150 °C	CHF 7'960.00
15.1880.16	Maxi 2 (2 chambers)	* 495 x 500 x 540 mm	* 156 l	-20 - +70 °C	CHF 6'001.00
15.1880.17	Maxi 2 Vitrum (2 chambers)	* 495 x 500 x 540 mm	* 156 l	-20 - +150 °C	CHF 7'357.00
15.1880.71	Mini-1 Desiccator Basic (1 compartment)	208 x 159 x 180 mm	6 l	-35 - +70 °C °C	CHF 252.00

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 10.05.2024 at 01:40 o'clock.