

# Cylindrical Tank with Cover and Spigot Nalgene

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Consumables](#) / [Containers](#) / Cylindrical Tank with Cover and Spigot Nalgene

## Product Data

LLDPE, translucent and lightweight

Economical tanks with graduations. With protruding rim (flange). Wall thickness ~2,4 mm.

- Good chemical resistance
- Heat resistant from -100 to +120 °C
- Conform to the requirements of USP VI and FDA 21 CFR 177.1520

Item No.	capacity	ext. D x height	Price
15.5410.54	19 l	280 x 380 mm	CHF 258.80
15.5410.55	28 l	300 x 460 mm	CHF 321.75
15.5410.56	38 l	330 x 510 mm	CHF 312.25
15.5410.57	57 l	360 x 710 mm	CHF 314.60
15.5410.58	113 l	460 x 760 mm	CHF 410.75
15.5410.59	280 l	560 x 910 mm	CHF 550.85

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 16.05.2024 at 02:45 o'clock.