

Erlenmeyer flasks with baffles DURAN

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Consumables](#) / [Containers](#) / [Flasks](#) / [Erlenmeyer flasks](#) / Erlenmeyer flasks with baffles DURAN

Product Data

Four baffles create turbulent flow, increasing the surface of the liquid and allowing greater gas exchange when compared to flat bottomed Erlenmeyer flasks.

- GL 45

Item No.	volume	diameter	Price
8.8336.01	250 ml	81mm	CHF 97.80
8.8344.01	500 ml	105mm	CHF 113.20
8.8354.01	1000 ml	136mm	CHF 38.25

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 30.05.2024 at 06:16 oclock.

www.huberlab.ch