

Clear Ampule, Pre-scored Wheaton

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Quality Assurance](#) / [Sample preparation](#) / [Glassware](#) / Clear Ampule, Pre-scored Wheaton

Product Data

- Pre-scored to eliminate filing
- Made from low extractable borosilicate glass that conforms to USP Type I and ASTM E 438 Type I, Class A requirements
- Stems can be pull or tip sealed

Item No.	vol.	d x h	Price
8.1767.72	1 ml	10,5 x 67 mm	CHF 84.38
8.1767.76	2 ml	12 x 75 mm	CHF 84.38
8.1767.79	5 ml	16,5 x 84 mm	CHF 130.32
8.1767.80	10 ml	19 x 107 mm	CHF 176.26
8.1767.82	20 ml	22,5 x 130 mm	CHF 259.78

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 17.05.2024 at 22:36 oclock.