

Serum Bottle, Graduated Wheaton DWK



Product in eShop

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Consumables](#) / [Containers](#) / [Flasks](#) / Serum Bottle, Graduated Wheaton DWK

Product Data

- Raised graduations indicate approximate capacities in milliliters
- Accepts standard stoppers and seals
- Manufactured from WHEATON 400 borosilicate glass that conforms to USP Type I requirement

Item No.	vol.	nech d	d x h	Price
8.2239.50	250 ml	15 x 30 mm	60 x 158 mm	CHF 358.75
8.2239.52	500 ml	15 x 30 mm	75 x 190 mm	CHF 325.20

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 03.05.2024 at 19:51 oclock.