

# Aluminum Seals, Open Top (Unlined) Wheaton

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Consumables](#) / [Containers](#) / [Flasks](#) / [Aluminum Seals, Open Top \(Unlined\) Wheaton](#)

## Product Data

- Wide range of styles and sizes
- Color selection includes natural, red, blue and green
- Open top seals expose stopper for sample retrieval with a syringe

Item No.	for	colour	Price
8.1562.24	Mouth OD: 11 mm	natural	CHF 179.30
8.1563.39	Mouth OD: 11 mm	green	CHF 205.15
8.1564.45	Mouth OD: 13 mm	natural	CHF 187.90
8.1565.51	Mouth OD: 13 mm	blue	CHF 216.25
8.1566.66	Mouth OD: 13 mm	red	CHF 216.25
8.1567.72	Mouth OD: 20 mm	natural	CHF 154.25
8.1568.87	Mouth OD: 20 mm	blue	CHF 205.15
8.1569.93	Mouth OD: 20 mm	red	CHF 205.15
8.1570.04	Mouth OD: 20 mm	green	CHF 205.15

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 17.05.2024 at 18:51 oclock.