

Headspace Vials, Crimp Top and flat Bottom Wheaton

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Chromatography](#) / [Vials & bottles](#) / [Headspace Vials, Crimp Top and flat Bottom Wheaton](#)

Product Data

- Manufactured from clear Type I borosilicate glass
- Accept 20 mm aluminum seals
- Variety of aluminum seal styles and septa materials available
- Purchase vials separately or together with seals in convenience packs

Item No.	size	description	d x h	Price
8.2885.59	6 ml	beveled finish	22 x 38 mm	CHF 64.25
8.2886.65	10 ml	beveled finish	23 x 46 mm	CHF 68.55
8.2887.71	20 ml	beveled finish	23 x 75 mm	CHF 87.05

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 17.05.2024 at 04:08 oclock.

www.huberlab.ch