

BOD Bottle Rack Wheaton



Product in eShop



[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Chromatography](#) / [Vials & bottles](#) / BOD Bottle Rack Wheaton

Product Data

- Conveniently and safely transports BOD bottles
- Racks are stackable
- Not autoclavable

Item No.	for	description	dimensions (l x w x h)	Price
8.2967.76	60 ml	Holds: 20 Bottles, max. Ø 49 mm	330 x 254 x 178 mm	CHF 152.90
8.2968.82	300 ml	Holds: 12 Bottles, max. Ø 76 mm	330 x 254 x 178 mm	CHF 124.35

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 17.05.2024 at 16:16 oclock.