

Samplepack YOUTILITY Labels DURAN® DWK



Product in eShop



Home / [Assortment](#) / [Quality Assurance](#) / [Sampling](#) / [Laboratory bottles](#) / Samplepack YOUTILITY Labels DURAN® DWK

Product Data

White, printable and self-adhesive

- Wide temperature performance range makes them ideal for use in fridges, freezers, autoclaves, incubators, and water baths
- Can be printed with office printers, copiers, or marked with technical lab pens
- The labels have a high adhesive strength, but are peelable and easy to remove which saves time during cleaning
- An free easy-to-use label creator web app for individual label design is available
- Chemically resistant to many typical laboratory chemicals, and solvents
- The labels are non-hazardous

Samplepack includes: [8.5827.94](#)

Item No.	Price
8.5817.99	for free

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 02.05.2024 at 20:45 oclock.