

## Sample Bags, Whirl-Pak® Speci-Sponge®



Product in eShop

[Home](#) / [Assortment](#) / [Medicine & Healthcare](#) / [Sampling](#) / [Sampling bags](#) / [Sample Bags, Whirl-Pak® Speci-Sponge®](#)

### Product Data

PE

Sterile, economical sampling bag with sponge is specifically designed for testing surfaces, equipment, and animal carcasses for the presence of *Listeria*, *Salmonella*, *E. coli*, and other food-borne pathogens. The 4x8 cm sponge comes in a Write-On Whirl-Pak bag. After adding diluent and wiping the surface with the moistened sponge, place it in the bag and whirl it closed to form a leakproof seal.

- Sterile, economical sampling bags

Item No.	capacity	thickness	l x w	Price
15.1129.18	1627 ml	0,102 mm	300 x 190 mm	CHF 350.90

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 18.05.2026 at 10:00 o'clock.