

Sapphire Low Profile 8-tube Strips 0.1 ml Greiner Bio-One

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Life Science](#) / [Genomics](#) / [Amplification](#) / [PCR strips](#) / [Sapphire Low Profile 8-tube Strips 0.1 ml Greiner Bio-One](#)

Product Data

Sapphire low profile PCR 8-tube strips with optically clear single attached flat caps for real-time PCR

Sapphire low profile PCR 8-tube strips with optically clear flat 8-cap strips for real-time PCR

- Optimal pour une PCR rapide
- Polypropylène de grade médical sans additifs interférents
- Fine paroi pour transfert de température homogène
- Épaisseur de paroi uniforme entraînant de faibles taux d'évaporation
- Adaptation universelle pour les thermocycleurs standard et RAPIDE
- Permet des volumes de réaction inférieurs à 20 l
- Efficacité, sans compromis de précision

Item No.	description	vol.	colour	Price
7.671 221	8-tube strips for qPCR, low profile, with individually attached optical flat caps	0,1 ml	natural	CHF 1'937.88
7.673 281	8-tube strips for qPCR, low profil, with optical flat caps	0,1 ml	natural	CHF 2'018.63

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 20.04.2026 at 11:10 o'clock.

www.huberlab.ch