


# CELLSTAR® TRI Suspension Culture Flasks Greiner Bio-One

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Life Science](#) / [Cell Biology](#) / [Cell culture consumables](#) / [Cell culture flasks](#) / CELLSTAR® TRI Suspension Culture Flasks Greiner Bio-One

## Product Data

### Standard Suspension Culture Flasks - Triple-Packed Products

- Triple-packed
- Hydrophobic surface, ideal for suspension cultures, hybridoma and embryonic stem cells
- Suspension culture flasks with standard screw cap (without filter)
- Canted neck
- Graduation on both sides

Item No.	Price
7.658 195-TRI	CHF 407.12
7.661 195-TRI	CHF 217.13

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 21.04.2026 at 06:36 o'clock.