


Accessories for Ion Exchange Stainless Steel Cartridges

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Water Purification](#) / [Ion exchanger & service](#) / Accessories for Ion Exchange Stainless Steel Cartridges

Product Data

Item No.	description	Price
16.1400.01	Pure water outlet, straight	CHF 273.80
16.1401.01	Pure water outlet, made of plastic, angled at 90°	CHF 273.80
16.1402.01	Distribution manifold: 1 inlet, 3 outlets	CHF 241.30
16.1402.02	Distribution manifold: 1 inlet, 2 outlets	CHF 175.00
16.1403.01	Shutoff tap	CHF 223.80
16.1690.01	2 x 2,5 Meter Schlauch R 3/4" max: 6 Bar	CHF 128.80
16.1690.06	---	CHF 63.65

More product images



* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 15.05.2024 at 19:36 oclock.

www.huberlab.ch