

## Media bottles tetragonal Greiner Bio-One

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Microbiology](#) / [Media preparation](#) / [Flasks, containers & tubes](#) / Media bottles tetragonal Greiner Bio-One

### Product Data

- Made of polyethylene terephthalate (PET)
- Radiation sterilized & with graduation
- Free of detectable DNase, RNase and human DNA
- Free of detectable endotoxins, free of cytotoxic substances
- Tetragonal
- Sterile

Item No.	vol.	Price
7.950 700	500 ml	CHF 132.85
7.950 710	500 ml	CHF 132.85
7.951 700	100 ml	CHF 163.47
7.951 710	100 ml	CHF 163.47
7.952 700	1000 ml	CHF 92.86
7.952 710	1000 ml	CHF 92.86

### More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 14.05.2024 at 02:42 oclock.

[www.huberlab.ch](http://www.huberlab.ch)