

## Cryo Container «Mister Frosty» Nalgene®



Product in eShop



Home / Assortment / Life Science / Microbiology / Media preparation / Cooling & freezing / Cryo Container «Mister Frosty» Nalgene®

### Product Data

Container made of PC, with blue HDPE lid, white HDPE tube holders and foam insert. Ensures the critical, repeatable cooling rate of 1 °C / minute required for successful cryopreservation and recovery.

**Easy to use**, requiring only isopropyl alcohol and a mechanical freezer

- Labeled with detailed step-by-step instructions.
- For 12 cryotubes or 18 cryotubes.
- Containers are stackable and can be stored at room temperature when filled with alcohol, so no freezer space is blocked
- Diameter: 117 mm
- Cap size: 120 mm

Item No.	for	Price
15.4177.06	18 Cryoröhrchen (1,0 - 2,0 ml)	CHF 168.60
15.5100.27	12 Cryoröhrchen (3,6 ml)	CHF 173.35
15.5100.28	12 Cryoröhrchen (4,5 - 5,0 ml)	CHF 173.35

### More product images



\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 01.05.2024 at 14:41 oclock.

[www.huberlab.ch](http://www.huberlab.ch)