

Micro hematocrit capillaries



Product in eShop



[Home](#) / [Assortment](#) / [Medicine & Healthcare](#) / [Hematology](#) / [Blood Collection](#) / [Micro capillary pipettes & pipetting aids](#) / [Micro hematocrit capillaries](#)

Product Data

Heparinised or non-heparinised, disposable, fire-polished ends; inside diameter equally over the total length (conicity 0.02 mm). Ungraduated, for a volume of 40 µl, for use in hematocrit centrifuges. Manufactured according to international specifications, with colour code. Heparinized over the entire internal surface, with sodium heparin, with a red colour code at one end.

- Packing: 100 pieces in a tube, 10 tubes in a box.

| Item No. | l x d | execution | Price |
|-------------|--------------------|---------------------------------|----------|
| 3.3530.04 | 75 x 1,3 / 1,4 mm | heparinised | CHF 5.28 |
| 3.3530.06 | 75 x 1,5/1,6 mm | heparinised | CHF 3.10 |
| 3.3530.08 | 100 x 1,5/1,6 mm | heparinised | CHF 9.90 |
| 3.3530.10 | 75 x 1,5/1,6 mm mm | heparinised, with mark at 60 mm | CHF 5.64 |
| 700.1602.13 | --- | --- | CHF 9.65 |

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 23.05.2024 at 20:33 oclock.