

PARAFILM® M sealing film

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Laboratory consumables & instruments](#) / [Consumables](#) / [Films](#) / [Parafilm](#) / PARAFILM® M sealing film

Product Data

stretches up to 200% and clings even around irregular shapes and surfaces. The sealing film is free of plasticizers and consists primarily of polyolefins and paraffin waxes

- It conforms to the provisions of the US FDA Food and Drug Administration when used below 55 °C and observance of GMP (Good Manufacturing Practice)
- Temperature range (continuous use): -45 °C to +50 °C
- Is resistant up to 48 hours against many polar substances, e.g. saline solutions, inorganic acids and alkaline solutions. After this period embrittlement may occur

Item No.	width	length	Price
15.1550.01	50 mm	75 m	CHF 23.20 CHF 58.00
15.1550.02	100 mm	38 m	CHF 17.75 CHF 39.00
15.1550.03	100 mm	75 m	CHF 46.50 CHF 104.00
15.1550.10	500 mm	15 m	CHF 48.00 CHF 120.00

Accessories

Item No.	Description	Price
15.1551.01	PARAFILM® M Cutter	from CHF 346.90
15.1551.03	PARAFILM DISPENSER	from CHF 11.95
15.1551.02	Acrylic dispenser for Parafilm sealing film	from CHF 205.65
15.1551.15	Box Top Parafilm® M Dispenser	from CHF 14.65

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 03.06.2025 at 07:18 oclock.

www.huberlab.ch