

## Beta safety box

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Safety & Hygiene](#) / [Disposal](#) / [Disposal containers](#) / [Beta safety box](#)

### Product Data

Made of acrylic, with tightly closing hinged lid, for contamination-free storage of vials

- Thickness 10 mm
- Without inserts

Item No.	description	Price
18.0003.02	inserts for 32 reaction vessels 1.5 + 2 ml	CHF 93.65
18.0003.03	inserts for 32 reaction vessels 0,5 ml	CHF 93.65
18.0003.04	Insert for 32 cryovials 2 ml	CHF 93.65

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 23.05.2026 at 01:33 o'clock.