


Cell Culture Dish CELLview with Glass Bottom for microscopy Greiner Bio-One

[Product in eShop](#)

[Home](#) / [Covid-19](#) / [Covid-19 Greiner Bio-One](#) / [GBO Research](#) / [Cell Culture Dish CELLview with Glass Bottom for microscopy Greiner Bio-One](#)

Product Data

The CELLview™ cell culture dish combines the convenience of a standard size 35 mm disposable plastic cell culture dish with the optical quality of glass, providing superior high-resolution microscopic images of in-vitro cultivated cultures.

It is made from high-grade polystyrene combined with an integrated glass bottom. The innovative design of the cell culture dish provides a single-plane, flat bottom with a consistent working distance and maximal planarity. Moreover, the dish bottom configuration facilitates optimal thermal conductivity and avoids thermal variations in heated platforms used for live cell imaging.

The subdivided version of the CELLview™ Dish enables simultaneous multiplex analyses of different cell lines, various stimulations or diverse transfections. Quartering the cell culture dish provides four individual compartments with a growth area of approximately 1.9 cm², allowing minimisation of cells and reagents required per individual assay.

Glass bottom features:

- High transparent achromatic borosilicate glass, hydrolytic class 1 (DIN ISO 719)
- Glass thickness 175 μm +/- 15 μm
- Maximal spectral transmission, no autofluorescence
- Exceptional planarity
- Manufactured according to ISO 8255-1:1986 (Optics and optical instruments – Microscopes – Cover glasses)

Advantages:

- Subdivided version enables simultaneous multiplex analysis
- Embedded glass bottom for maximal planarity
- Additional TC surface treatment and Advanced TC™ surface modification available

Applications:

- Phase contrast microscopy
- Fluorescence microscopy
- Confocal microscopy
- Live cell imaging
- Differential interference contrast microscopy
- Polarised light microscopy
- Fluorescence-in-situ hybridisation (FISH)

*) 0.1 ml for seeding or staining only on glass area, 0.5 ml for cultivation in the complete compartment

Item No.	description	d x h	area	working volume	TC treated	Price
7.627 860	1 compartment	35 x 10 mm	8,7 cm ²	5,0 ml	+	CHF 268.00
7.627 870	4 compartments	35 x 10 mm	1,9 cm ²	0,1 / 0,5* ml	+	CHF 544.00
7.627 871	4 compartments	35 x 10 mm	1,9 cm ²	0,1 / 0,5* ml	-	CHF 265.00 CHF 530.00
7.627 965	1 compartment	35 x 10 mm	8,7 cm ²	5,0 ml	Advanced TC	CHF 568.00
7.627 975	4 compartments	35 x 10 mm	1,9 cm ²	0,1 / 0,5* ml	Advanced TC	CHF 1'140.00

* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

** Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

*** Further information such as technical information and safety data sheets can be found online in our eShop.

**** The PDF file was created on www.huberlab.ch on 21.04.2026 at 06:36 o'clock.