

## Pipette baskets

[Product in eShop](#)

[Home](#) / [Assortment](#) / [Research & Development](#) / [Liquid Handling](#) / [Pipetting](#) / [Glass pipettes](#) / [Pipette baskets](#)

### Product Data

Out of PE-HD, for pipette rinsing systems.

Pipette rinsing systems:

Trouble-free and gentle cleaning of pipettes with repetitive siphon rinsing. Water supply connection via hose nozzle, outer-Ø 12 mm. The complete cleaning system consists of a rinsing apparatus (15.1556.00 & 15.1556.99), soaking jar (15.1556.10 & 15.1556.12) and pipette basket. Please order each component separat

Item No.	for pipette length	description	foot Ø	height	Price
15.1556.15	460 mm	Basket height incl. handle: 645 mm	145 mm	280 mm	CHF 82.45
15.1556.16	600 mm	Basket height incl. handle: 870 mm	145 mm	280 mm	CHF 79.95

\* The prices are non-binding and are to be understood as selling prices in Swiss francs without value added tax (VAT), as well as all other fees, charges and taxes. The prices displayed in the eShop may differ from the PDF file due to regular updates.

\*\* Please note that when ordering chemicals and detergents, transport and packaging costs for hazardous goods as well as legally prescribed fees are charged. These will be shown in detail on the order confirmation, which you will receive in addition to the confirmation of receipt.

\*\*\* Further information such as technical information and safety data sheets can be found online in our eShop.

\*\*\*\* The PDF file was created on [www.huberlab.ch](http://www.huberlab.ch) on 11.05.2024 at 17:23 oclock.